miniland

nebuClean

Español	06	Português	31	Deutsch	59	Polski	87
English	19	Français	45	Italiano	73	Ελληνικά	100

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INDEX

1. Introduction
2. Contents
3. Important safety precautions
4. Product features
5. Instructions
6. Maintenance
7. Technical specifications
8. Product disposal information

1. INTRODUCTION

We congratulate you on acquiring this portable mesh nebulizer.

Before using this device, please carefully read all instructions in this manual. Keep the manual for future reference as it contains important information.

The features described in this user guide are subject to change without notice.

2. CONTENTS

- 1 main unit
- 1 medication cup (piezoelectric component)
- 1 infant mask
- 1 adult mask
- 1 mouth piece
- 1 USB cable
- 1 carrying bag
- 1 instruction manual and warranty

If any of the above-mentioned contents are missing, please contact your point of sale.

3. IMPORTANT SAFETY PRECAUTIONS

WARNING

- Use this appliance under the guidance of a professional doctor. For the type, dose, and regimen of the medicine, be sure to follow the instructions of a doctor.
- Children or disabled patients should be assisted by adults or medical personnel when using this appliance.
- Keep this appliance away from the reach of children.
- Do not use this appliance under flammable gases.
- Use this appliance only for its intended use as described in this manual, and do not use any

attachments or accessories not recommended by the manufacturer.

- Cut off the power when the appliance is abnormal.
- No modification of this equipment is allowed.
- This equipment is not suitable for high suspended liquids or viscous liquids.
- This equipment is intended for use only in conscious patients.
- During the use of this appliance, if any irregularity or discomfort is observed, stop and seek medical attention immediately.
- This appliance is not intended for life support.
- Do not use this appliance near high frequency, microwave ovens and mobile phones.
- Prior to the initial operation, or if the appliance has not been used for a long time, proceed with the cleaning and sterilization procedures described in this user manual.

PRECAUTIONS

- Use and store the appliance only in the operating and storage environments specified in this manual.
- Do not use this appliance if the medication cup doesn't contain any medicine.
- Different patients should use individual clean medication cups and mouth pieces to avoid cross infection.
- Do not pour water inside the device.
- The main unit, medication cup, and nebulizer accessories should not be placed in boiling water for disinfection and should not be placed in a microwave oven for drying.
- The product has an in-built lithium battery, there is no need to replace it.
- Please pay attention to this device when charging, do not overcharge.

CONTRAINDICATIONS

- The patients with serious lack of oxygen or respiratory failure are prohibited to use this appliance.

IMPORTANT: In case of product quality issues, contact Miniland's customer service. Do not disassemble or repair this product yourself. Changes and modifications to the appliance not made by an authorized service center are not covered by the product warranty.

4. PRODUCT FEATURES

4.1. WORKING PRINCIPLE

This product is mainly composed of piezoelectric components. The piezoelectric components convert electrical energy into mechanical energy and generate ultrasonic vibration. The shock wave squeezes the liquid in the medicine cup, so that the liquid medicine passes through the micropores on the atomizing sheet to spray the liquid medicine. It is sprayed out as mist and used by inhalers.

4.2. INTENDED FOR USE

Atomize the medicine, for treatment. The product is used for the conscious patients. It can be used on all children and adults.

4.3. DESCRIPTION OF PARTS (FIGURES 1, 2 AND 3)

1. Cup lid
2. Medication cup
3. Mesh
4. Cup buckle
5. Medication cup button
6. USB plug (Type C)
7. Electrode
8. Medication cup indicator
9. Working indicator
10. Low voltage indicator
11. Power button
12. Medicación cup protection cover
13. Mouth piece
14. Mask

5. INSTRUCTIONS

5.1. INSTRUCTIONS FOR USE

Make sure the nebulizer components are clean, hygienic and dry before use. After, follow the next steps:

1. Add the liquid medicine or saline or hypertonic serum (figures 4 and 5):

- a. Push the medication cup to the front of the main unit after pressing the back cup button (5) on the main unit.
- b. Open the medication cup cover and pour the medicine into the medication cup and close the cover completely:

- Add at least 1 ml of liquid, and do not add medicine over the indicated line marked on the medication cup (2).

- Remove the medication cup (2) from the main unit when adding medicine.

- Fasten the medication cup lid to prevent the liquid from leaking out.

- Do not shake the nebulizer violently when the medication cup is filled with liquid.

c. Put the cup back into the main unit. You will hear a "click" sound once the cup is in the correct position.

Please keep the electrode on the main unit and the medication cup clean (no dirt or liquid residue) and confirm that the medication cup is assembled in place. Otherwise, the electrode may be abnormally turned on and the nebulizer may not work properly.

2. Install the mask or mouthpiece (figure 6):

- a. Select the appropriate accessory (mouthpiece, adult mask or child mask) according to your needs. Check if it is clean, hygienic and dry before first use.
- b. Plug the accessory into the medication cup.

Note 1: Please, use the inhalation device according to the doctor's instructions.

Note 2: Apart from the accessories and parts that are sold separately, use masks and mouthpiece recommended by the manufacturer only.

3. Turn on the device pressing the power button (11). The appliance will start working and the work indicator (9) will light in green.

After the power is turned on, the nebulizer will undergo a short start-up period (within about 2 seconds) and then spray normally.

During the spray process, there may be a noise caused by high frequency vibration. Please shake the nebulizer gently to reduce the probability of sound.

If there is no medicine or very poorly conductive liquid (for example, distilled water) inside the medication cup, the medication cup indicator will flash in blue, and the device will turn off automatically.

Cleaning mode

In the shutdown state, if the user keeps the power button (11) pressed for 3s, the "cleaning mode" will activate. In this mode, the green work indicator will light slowly on and off. The device will work for 3 minutes and then automatically shut down.

4. Hold the nebulizer smoothly, relax and start breathing

- Since the characteristics of each liquid are different, when the nebulizer can't automatically shut down after some of the liquid medicine is exhausted, press the power button to shut down the appliance immediately to avoid damage to the mesh.
- When the liquid is about to be exhausted, it is recommended to tilt the front side of the nebulizer (on the side of the spray direction) slightly towards the patient, which is good for contacting the remaining liquid to the atomized piece for nebulization.
- Do not cover the vents on the mask and mouthpiece with your hands or other objects during the atomization process.
- During the use of this appliance, if any irregularity or discomfort is observed, please stop the use and seek medical attention immediately.
- During the atomization process of some liquids, a large amount of foam will be concentrated in the vicinity of mesh in the medication cup, which may easily cause the spray to be interrupted. In this case, turn off the nebulizer immediately, shake the nebulizer slightly, and then restart.

5. Power off the device

During use, if you need to stop the device, you can press the power button to turn it off and the work indicator (green light) will turn off. Otherwise, the nebulizer will turn off after 20 minutes of continuous operation or after the liquid medicine is exhausted.

In this case, the nebulizer will emit a high-frequency sound, the medication cup indicator (blue light) will flash, and the nebulizer shut down automatically.

5.2. BATTERY CHARGING

This device comes with an inbuilt lithium battery. It's recommended to connect it to the power supply to activate the battery for first time use or after a long period of non-use.

To charge the battery, open the Type-C USB socket protection plug (6) on the back of the main unit, insert the Type-C USB plug into the socket hole, and then plug the adapter power plug into the electrical outlet.

The work indicator will start flashing in green indicating that the charging has begun. Once is completed, the light will stop flashing.

After charging is completed, unplug the adapter power plug and the Type-C USB plug, and close the Type-C socket protection plug (6).

If the battery is about to run out, the low voltage indicator (10) will flash in yellow. At this point, follow the previous indications to charge it.

Remember to unplug the device from the electrical outlet after charging is completed, do not connect the power supply for long time.

Important remarks:

- Do not use the nebulizer when charging.
- It is recommended to charge the lithium battery once when the product is not used for three months.
- It is recommended to use the specified power adapter for this product.

5.3. TROUBLESHOOTING

FAILURE	CAUSE OF FAILURE	COUNTERMEASURES
The device doesn't turn on	Low battery	Charge the lithium battery and reboot after the charge is completed.
	The electrode on the main unit or the cup is dirty	Clean the electrode and reboot.
The device is on but there is no operation	If the low voltage indicator (10) is flashing, the device doesn't have enough battery	Charge the lithium battery and reboot after completed charging.
	The mesh (3) is dirty or blocked	Clean the medication cup according to the instructions. Replace the medication cup if the fault is not excluded after cleaning.
	Medication cup is not assembled in place	Reassemble the medication cup and reboot.
	The medication cup has been filled with liquid medicine for a long time and a water film has formed	Clean the medication cup and reboot.
The device turns on, the medication cup indicator (8) flashes and the device turns off	There is not enough liquid on the medication cup	Add medicine or saline.
	Very poorly conductive liquid inside medication cup (such as distilled water)	Add medicine or saline.
	Medication cup is not assembled in place	Reassemble medication cup, and reboot.
	The electrode on the main unit or the cup is dirty	Clean the electrode and reboot.
	The liquid has not been in contact with the mesh for more than about 8 seconds	Tilt the front of the nebulizer to the user to make the liquid contact the mesh.
Auto-off during use	Liquid is exhausted	Add medicine or saline.
	The liquid is not in full contact with the mesh	Tilt the front of the nebulizer to the user to make the liquid contact the mesh.
	There is a liquid residue on the mesh, and a very weak conductive liquid (such as distilled water) is placed in the medicine cup	Add medicine or saline.

Liquid is exhausted, but the nebulizer doesn't turn off	The conductivity of the liquid is too strong due to its chemical properties	Manually shutdown the device.
	Some liquids may produce some foam in the medication cup	Clean up the foam and reboot.
Liquid leaks	The liquid on the medication cup exceeds the maximum volume	Pour out excess liquid and restart.
	During use, the medication cup is being shaken	Hold the nebulizer steadily during use.
	The medication cup is damaged, or the sealed silicone ring has aged	Replace the medication cup for a new one.

If the nebulizer still does not work properly after taking the above measures, we recommend that you contact Miniland.

6. MAINTENANCE

6.1. CLEAN

6.1.1. CLEANING RESIDUAL LIQUID FROM THE MEDICINE CUP

- Open the medication cup and discard all the remaining liquid.
- Pour in some clean water and close the medication cup completely. Gently shake the cup to dissolve the residue in purified water.
- Pour off the liquid and refill with some clean water.
- Long press the power button key for about 3 seconds to enter the "cleaning mode". The device will work for 3 minutes and then automatically shut down to clean the mesh.
- Put on the protective cover on the medication cup.

Notes:

- Ensure to clean the liquid completely after each use, otherwise the mesh could block.
- After entering the cleaning mode, if the nebulizer is running out of water and there is no spray, please turn it off immediately.
- If the residual liquid has not been removed, it can be washed with warm water.
- Clean water includes distilled water and pure drinking water.

6.1.2. CLEANING THE MAIN UNIT AND ACCESSORIES

Accessories:

- Remove the medication cup and mask or mouthpiece.
- Rinse the cup, mask or mouthpiece with clean water. Or soak in clean water for 15 minutes.

Main unit (figure 7):

- Clean the surface of the main unit with a piece of wet fine cotton fabric and dry it with a dry cotton fabric.
- Clean the electrodes on the main unit and the cup assembly to ensure that the cup is securely attached to the main unit.
- Wipe the cleaned parts with a clean gauze and dry thoroughly.
- Install the medication cup on the main unit, and store in a clean place.
 - Do not wipe the nebulizer with volatile corrosive liquids such as benzene, benzine or thinner.
 - Unplug the power adapter from the socket after turning off the power.
 - To prevent residual paper dust or cloth from entering the cup and clogging the atomizer, do not use cotton paper or easy-to-hair cloth to wipe the cup assembly.
 - Do not clean the components of the product in the dishwasher, and do not dry the product components in the microwave oven.

6.2. STERILIZE

Please disinfect the mouthpiece and mask after each use. If the parts are seriously

polluted, please replace the parts in time.
It is recommended to follow these disinfection methods:

- a) Alcohol disinfection: disinfect the nozzle and mask with 75% medical alcohol.
- b) After disinfection, rinse thoroughly with clean water to avoid residual disinfectant.
- Do not disinfect with a sub-chlorinated disinfectant. Such disinfectant may cause damage to the components.
- If you use disinfectant and it accidentally gets into your eyes or sticks to the skin, rinse

with clean water as soon as possible or go to a nearby medical facility.

6.3. STORAGE

After each use and cleaning and disinfection, confirm that the body and the cup are dry and store them in the transportation bag. Place the bag in a dry and cool place and keep it away from the reach of children.

7. TECHNICAL SPECIFICATIONS

7.1. DEFINITION OF SYMBOLS


	Fragile, handle with care
	Keep dry
	This way up
	Stacking layer
	Do not re-use
	Medical device
	Caution
	Type BF applied part
	Manufacturer
	Authorized representative in the European Community
IP22	IP code of the device: this device's grade of against ingress of solid foreign objects -- $\geq 12.5\text{mm}$ diameter (and the against access to hazardous parts with finger); the grade of waterproof is dripping (15° tilted)
	Complies with the European Medical Device Directive 93/42/EEC and amended Directive 2007/47/EC. Notified Body is SGS

	Date of manufacture
	Batch code
	Follow operating instructions The safety colour of this safety sign is blue, and the graphical symbol colour is white.
	Disposal in accordance with Directive 2012/19/EU (WEEE)

7.2. TECHNICAL INFORMATION

Device name	Portable (ultrasonic) nebulizer
Model	89651 nebuClean (NBM-6B)
Power supply	d.c.3.7V lithium battery
Power consumption	<2VA
Standby power consumption	<0.02VA
Nebulization rate	>0.2mL/min
Particle Size (MMAD)	<5µm
Amount of residue	≤0.5ml
Noise	≤50dB(A)
Mesh oscillation frequency	110 ± 10KHz
Battery life	>2 hours
Time AUTO-OFF	About 20 minutes
Operation environment	Temperature range: 10~40°C Humidity range: 10%-95%RH (Non-condensing) Atmospheric pressure range: 860-1060hPa
Storage/Transportation environment	Temperature range: -20~70°C Humidity range: 10%-95%RH (Non-condensing) Atmospheric pressure range: 500-1060hPa
Type of protection against electric shock	Internal power supply Type BF applied part
Grade of waterproof	IP22
Service life of main unit	2 years
Dimensions	Approx 42(L)*63(W)*113(H)mm
Weight	Approx 100g
Last update of the manual	October 2024

The equivalent volume particle size distribution curve of the atomized particles is as follows:

Test conditions: ambient temperature $25 \pm ^\circ\text{C}$, relative humidity $\leq 85\%$, using pure water test. (Test data may vary depending on test conditions and changes in liquid medicine)

7.3. EMC INFORMATION

The device needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and the device can be affected by portable and mobile RF communications equipment.

EMC Warning:

- 1) The portable (ultrasonic) nebulizer is intended to be used in the hospital and family setting except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields near the unit. This may result in incorrect operation of the unit.
- 3) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation.
- 4) Caution: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.

- 5) Use of accessories other than those specified or provided by the manufacturer of this the nebulizer could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Table 1

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The device is suitable for use in all establishments including those directly connected to the public low-voltage power supply network that supplies power to buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations and flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and Manufacturer's Declaration — Electromagnetic Immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle g) At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycle	0% UT; 0,5 cycle g) At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Table 3

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz	3 Vrms 150 kHz to 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 385MHz-5785M Hz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment	10 V/m 80 MHz to 2.7 GHz 385MHz-5785 MHz Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment	Recommended separation distance $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.7 GHz where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,a should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol: 
---------------------------------	---	---	---

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a) The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

c) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device				
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output of transmitter W	Separation distance according to frequency of transmitter			
	150 kHz to 80 MHz outside ISM and amateur radio bands	150 kHz to 80 MHz in ISM and amateur radio bands	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{12}{V_2} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.20	0.035	0.07
0.1	0.38	0.63	0.11	0.22
1	1.2	2.00	0.35	0.70
10	3.8	6.32	1.10	2.21
100	12	20.00	35	70
For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

8. PRODUCT DISPOSAL INFORMATION

- Remember to remove the batteries and dispose of them in a responsible manner when disposing of the baby monitor.
- You should not mix products bearing the crossed-out rubbish bin symbol with your domestic rubbish.
- These products should be taken to the collection points allocated by the local authorities for appropriate treatment.

- Alternatively, you may contact the point of sale where you purchased the product.

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© Miniland, S. A. 2025
Dongguan SIMZO Electronic Technology Co., Ltd.
MedPath GmbH