
miniland

sweetBeat

Español	13	Português	23	Deutsch	33	Polski	44
English	18	Français	28	Italiano	39	Русский	49

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TABLE OF CONTENTS

1. Introduction
2. Content
3. Safety instructions
4. Product characteristics
5. Instructions for use
6. Maintenance
7. Technical specifications
8. Identifying and solving problems
9. Information on disposing of the product

■ 1. INTRODUCTION

We would like to congratulate you on having acquired our new innovative Doppler ultrasound heartbeat monitor which will enable you to listen to your baby's heartbeat as of 16 weeks of pregnancy on your mobile device using Bluetooth technology. Read the instructions set forth in this manual carefully prior to using the product. Keep it for use as a reference in the future as it contains important information.

The characteristics described in this manual are subject to change with no prior warning.

■ 2. CONTENT

- 1 heartbeat monitor
- 1 USB cable charger
- 1 instructions manual and guarantee

If any of the aforementioned contents are missing, please contact your point of sale.

■ 3. SAFETY PRECAUTIONS

This device is rechargeable, energy-independent and complies with EU safety standards in relation to healthcare products. The following specifications should be complied with in order to avoid any potential damage.

1. The device is not suitable for use in environments containing inflammable and/or anesthetic gas as this could cause explosions.
2. Do not throw the device on to a fire in order to prevent it from exploding.
3. Do not use the heartbeat monitor on more than one patient simultaneously.
4. Do not use this device if there is high-frequency surgical or magnetic resonance equipment in use. This could harm both the operator and the patient. The correct operating distance between the device and the mobile player is 5 m.
5. This device has been designed to operate on a continuous basis. Although the device is manufactured splash-proof, avoid contact with liquids. Keep the device clean and keep it away from vibrations.
6. Do not sterilize the equipments in high-temperature environments. Also avoid the use of electronic, ultraviolet rays, etc.
7. Prior to use, check for any damage to the equipment which might place the safety of the user at risk. This check should be

conducted at least once a month to ensure the device can be repaired prior to use in the event a problem is detected.

8. If there is doubt as to fetal well-being after using the unit, consult with your doctor.
9. This device is a portable fetal heartbeat detector. Under no circumstances whatsoever may it replace medical fetal monitoring. We recommend the end user seek basic professional healthcare guidance prior to using the device for the first time.
10. We recommended the minimum use of this equipment immediately after a prescribed clinical monitoring procedure.

■ 4. CHARACTERISTICS OF THE PRODUCT

4.1. DIAGRAM OF THE PARTS OF THE PRODUCT (FIGURE A)

- 1) Power and fetal heartbeat indicator
- 2) 3.5 mm audio jack port
- 3) Mini-USB charge port
- 4) On / Off button
- 5) Transducer (lower part)

■ 5. INSTRUCTIONS FOR USE

5.1. CHARGING THE BATTERY

This heartbeat monitor is equipped with a rechargeable battery. When the battery is flat the power indicator (1) will light up orange. In order to recharge the battery:

- 1) Connect the end of the charge cable micro-USB supplied to the micro-USB port located on the side of the device covered with a small piece of silicone (3). Connect the other end to your computer. The device will then begin to charge.
- 2) The power indicator will flash orange while the device is charging. Once charged it will flash green.

Note 1: Do not attempt to remove the battery from the apparatus.

Note 2: You may not use the apparatus to read a fetal heartbeat while it is charging.

5.2. DETECTING THE HEARTBEAT

As a means of ensuring your experience with the heartbeat monitor is as satisfactory as possible, the transducer must be placed in the correct position on the abdomen. As such, feel for the first third of the vertical area joining the navel and the pubis and apply a generous amount of conductive gel.

Then, depending on whether the patient is in phase 1 (16-24 weeks), phase 2 (24-32 weeks) or phase 3 (32-40 weeks) of pregnancy, feel for the heartbeat of the fetus as illustrated in **Figure B**.

Obviously, this area will expand over the pregnancy. This is due to the fact that as the fetus

grows, it, and consequently the heart, moves upwards and the mobility towards the sides increases. The ideal position for clearly detecting the fetal heartbeat is when the apparatus is pointing towards the rear of the fetus' left shoulder. Ensure the transducer (5) is always fully in contact with the skin.

Once the best position for listening to the baby's heartbeat has been found, you can use your mobile device with the eMyBaby application to monitor the fetal heartbeat. Or if you prefer, you can also listen using the headphone connection in the 3.5 mm jack input (2).

If the fetus moves position during the examination, you should change the direction of your heartbeat monitor until you can hear the cardiac activity with the greatest possible clarity.

Note 1: Miniland has a special gel to favor the conduction of ultrasound, whose reference is 89211. You can purchase it through our online store (<https://www.minilandgroup.com/>). Also, this type of gel can be found in pharmacies and parapharmacies.

Note 2: You should take into account that the power of the ultrasound signal of this device is far lower than that of standard industrial devices, meaning it is much safer for the fetus. However, its power of detection is less efficient for the same reason and despite the high sensitivity of the transducer. Hence, you need to find the best position for detecting the fetal heartbeat in order to obtain the greatest possible precision.

Note 3: It is not always possible to hear a baby's heartbeat before 16 weeks and will depend largely on various factors. Among them, the height and figure of the mother. Mothers who are taller or overweight may find it more difficult to hear their baby's heartbeats for the first few weeks.

5.3. USE OF THE HEARTBEAT MONITOR

Please follow the steps listed below in order to connect your mobile device to sweetBeat and make use of the same:

1. Turn on the heartbeat monitor by pressing the on/off button (4) for three seconds until the power indicator (1) lights up green.
2. Download and install eMyBaby in accordance with the instructions provided at the beginning of this manual.
3. Access the eMyBaby "My Pregnancy" section and the "Heart Listener" sub-section (Figure 1). Once you have accessed this section the mobile device on which you have installed the application will start searching for bluetooth on your sweetBeat.

Note 1: Ensure the Bluetooth system on your mobile device is turned on prior to starting the process.

4. The icon of your heartbeat monitor device should appear on the screen to enable you to complete the connection by selecting it. Once this process has been completed, the device will appear on the screen as connected (Figure 2). This process will be repeated each time you wish to use the apparatus.
5. Place a generous amount of gel in the transducer.
6. On starting to use the device, you should find the best position to point the transducer towards the fetus' heart in accordance with the instructions provided in section 5.2. On placing the device on the mother's belly it will emit a green light while the fetal heartbeat registers between 110 and 160 beats per minute, and an orange light when beyond this range.
7. Once the sound begins to come through clearly, a graph of the baby's heartbeat will appear in real time (Figure 3). If you wish, you can also register the baby's movements manually while reading its heartbeat.

Note 2: You may turn up or turn down the volume of the sounds detected using the control on your mobile device or by using the corresponding icon on the screen.

8. If you wish, you can record the sounds and the graph by pressing the "Record" button. You may access these files at a later date in the "Your Sounds" sub-section.

Note 3: The device turns itself off automatically when it has not received a signal for four minutes as a means of increasing the life of the battery.

5.4. ADDITIONAL SECTIONS OF THE "My Pregnancy" SECTION ON eMyBaby (FIGURE 1)

In addition to the possibility of using the heartbeat monitor with eMyBaby as described in the previous section, the application features many more functions to enable you to follow the entire pregnancy from close up and in detail. The following additional sections are specifically geared to enabling parents to monitor the entire process.

5.4.1. Your sounds

This section registers the sounds recorded on your heartbeat monitor (explained in item 5.3 of this manual), indicating the date and time of the recording. Moreover, you can share these sounds with your loved ones.

ENGLISH

5.4.2. Baby

This section is specifically intended to register data in relation to the growth of your baby.

You can choose from the following:

- Insert your baby's ultrasounds (Figure 5) with comments.
- Insert the data provided by your doctor on different measurements of the fetus (Figure 6), which are then summarized in chronological order.
- Obtain graphs on this data showing your baby's progress. (Figure 7).

5.4.3. Mum

This section enables you to register important values in relation to the changes the mum undergoes during the pregnancy.

- Evolution of weight: In this section the mother can register her weight month by month, obtaining a graph showing progress (Figure 8).
- Blood pressure: To record each reading and to monitor the same (Figure 9).
- Fetal movements: This section shows if the baby has moved or not during the gynecological monitoring. (Figure 10).
- Height of the womb: This data is normally monitored in detail. The application allows you to register the values taken by your gynecologist and summarizes the progress (Figure 11).
- Stomach: A good test of the evolution of your future baby is the growth of the stomach. You can register photos of the same with comments on each one, and then edit or delete them (Figure 12).

5.4.4. Doctor's appointments

The application also enables you to keep a record of your doctor's appointments divided into specialist areas (Figure 13). You can add them selecting the respective specialist area and entering the data on the same, which will then appear in the "Calendar" section on the eMyBaby application.

5.4.5. Contractions

This section (Figure 14), at the time of giving birth, will enable you to record each contraction and the duration of each. This information can be summarized as illustrated in figure 14.

5.4.6. Baby's diary.

With the aim of enabling you to register all the events of this important phase for parents and the baby, the application features a diary in which you can register all the events occurring during the pregnancy, and even stores multimedia files (Figure 15).

6. MAINTENANCE

- This device is a precision tool and needs to be handled with care. All remains of gel from previous uses should be removed from the contact area prior to use. This will help improve the results of the exam and increase the service life of the equipment.
- Prior to using the heartbeat monitor, ensure there is no obvious damage which might affect its performance, as this could also cause harm to the user. If the equipment is damaged dispose of it and replace it.
- The accuracy of the measurement of the fetal heartbeat cannot be adjusted manually. In the event of doubt with regard to this accuracy, the user may compare the same by using other medical devices such as a stethoscope.

6.1. Cleaning

- Keep the equipment clean and free of dust by ensuring regular cleaning with a clean smooth and dry cloth.
- If necessary, you can also use a cloth dampened with water or soapy water.
- Do not use aggressive solvents such as acetone for cleaning.

ATTENTION: The use of abrasive materials such as steel wool and similar is prohibited.

Do not submerge the apparatus in any liquid substance.

Do not leave any cleaning fluid on the surface of the apparatus.

Note: The probe should be cleaned with a 70% alcohol solution or isopropyl alcohol. Then leave it to dry or use a cloth.

6.2. Disinfection

As explained previously, the top of the equipment and the transducer should be cleaned with a 70% alcohol solution or isopropyl alcohol, and then dried with a cloth to remove any debris from the apparatus. No other kind of cleaning or disinfection is allowed in order to prevent the device from being damaged or malfunctioning.

6.3. Transportation and storage

This device should not be stored in erosive, moist or poorly ventilated environments or in extreme temperatures. Moreover, the equipment should be protected against impact or powerful vibrations when being transported.

7. TECHNICAL SPECIFICATIONS

- Bluetooth standard: 4.0
- Compatible with Bluetooth standard: 3.0 and 4.0

7.1. Physical Characteristics

- Size: 80 mm (diameter) x 43 mm (height)
- Weight: 75 g ± 10 g (including the battery)

7.2. Operating environment

- Temperature: From +5 °C to +40 °C
- Humidity: ≤ 80%
- Atmospheric pressure: From 86 kPa to 106 kPa
- Service life: 5 years

7.3. Transportation and storage

- Temperature: From +10 °C to +55 °C
- Humidity: <93%
- Atmospheric pressure: From 86 kPa to 106 kPa

7.4. Measurement of the fetal heartbeat

- Sensitivity: ≥ 90 dB
- Measuring range: From 50 to 210 beats per minute
- Accuracy: ± 2 beats per minute

7.5. Electrical specifications

- Battery: Rechargeable 3.7 V Ion-Lithium

7.6. Specifications of the probe

- Type: Continuous Doppler
- Water resistance: IPX1
- Acoustic operating frequency: 2.5 MHz ± 15%
- Ultrasonic output power: ≤ 20 mW
- Effective emission area: 245 mm² ± 15%
- P: < 1 MPa
- Palob: 20 mW/cm²
- Ispta: < 100 mW/cm²

NORMALIZED SYMBOLS

- The name and the address of the Manufacturer

- Authorized Representative in the European Community

- Batch code

- Attention: Read the operating Instructions for use!

- Applied part of type B

- Disposal in accordance with Directive 2012/19/UE (WEEE)

- Complies with the European Medical Device Directive (93/42/EEC) and amended by directive 2007/47/EC requirements.

0482

Notified body: MEDCERT GmbH (CE0482)

- Serial number

8. IDENTIFYING AND SOLVING PROBLEMS

Problem	Possible cause	Solution
Nothing happens when you try to turn the apparatus on	Flat battery	Charge the heartbeat monitor as illustrated in item 5.1.
You mobile device is unable to detect the heartbeat monitor in its range of coverage when conducting the Bluetooth device search	The heartbeat monitor is not turned on	Turn it on as illustrated in item 5.3.
	Your mobile Bluetooth system is not turned on	Turn on your mobile Bluetooth system as illustrated in item 5.3
	The heartbeat monitor battery has gone flat	Charge the heartbeat monitor as illustrated in item 5.1.
	You mobile device is not compatible with Bluetooth BLE 4.0	Check sweetBeat compatibility in section 7. TECHNICAL SPECIFICATIONS
	sweetBeat is out of the range of your mobile device's Bluetooth	Move your mobile device closer to the heartbeat monitor

ENGLISH

Problem	Possible cause	Solution
The removable module does not play any sound	The volume on your device is too low	Turn up the volume on the device
	The heartbeat monitor is not turned on	Turn it on as illustrated in item 5.3.
	The heartbeat monitor battery has gone flat	Charge the heartbeat monitor as illustrated in item 5.1.
The reading is not exact or you have doubts in relation to the result	Check the transducer is clean	Clean the transducer as explained in section 6.1 of this manual
	Check the means of measurement is correct	Ensure you have read the manual and know how to use the heartbeat monitor correctly
Other faults	Contact the MINILAND technical assistance department	

9. PRODUCT DISPOSAL INFORMATION

- Remember to dispose of the product responsibly.
- Do not mix products bearing the crossed -out bin symbol with your general household waste.
- For the correct collection and treatment of these products take them to the collection points designated by your local authority.
Alternatively contact the retailer who sold you the product.

Hereby, Miniland S.A. declares that this "89304 sweetBeat" is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

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© Miniland, S.A. 2019
Shenzhen Luckcome Technology Inc., Ltd.