
miniland

thermosense plus


Español	5	Français	41	Polski	77
English	17	Deutsch	53	Ελληνικά	89
Português	29	Italiano	65	Русский	101


1 Figura • Figure • Abbildung


2 Figura • Figure • Abbildung


3 Figura • Figure • Abbildung


4 Figura • Figure
Abbildung


5 Figura • Figure
Abbildung

6 Figura • Figure
Abbildung

7 Figura • Figure
Abbildung

ÍNDICE

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INDEX

1. Introduction
2. Contents
3. Warnings and precautions
4. Product information
5. Instructions for use
6. Product care
7. Technical specifications
8. Information on the disposal of the batteries and the product

1. INTRODUCTION

Thank you for purchasing this forehead & ear infrared thermometer. This thermometer can be used to read the body temperature by measuring the ear and forehead temperature, which is suitable for adult and children (ear test mode only for the child above 3 months).

Before using this device, carefully read all the instructions given in this manual. For safe and correct use of this product, be sure to read and fully understand the safety precautions explained in this user manual. Please keep in a safe place as it contains important information.

The features described in this user's guide are subject to modifications without prior notice.

2. CONTENTS

- 1 Infrared thermometer
- 1 Pouch
- 2 AAA Batteries
- 1 Operating manual and guarantee card

If any of the above is missing, please contact your retailer.

3. WARNINGS AND PRECAUTIONS

1. Keep out of reach of children under 12 years.
2. Never immerse the thermometer into water or other liquids (not waterproof). For cleaning and disinfecting, please follow the instructions in the "Care and cleaning" section.
3. Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using it on children.
4. Keep the thermometer away from direct exposure to the sun and keep it in a dust-free, dry area, well-ventilated place at a temperature between 5°C (41 °F) - 40°C (104°F). Do not use the thermometer in high humidity environments. (>95% RH)
5. Do not use the thermometer if there are signs of damage on the measuring sensor or on the instrument itself. If damaged, do not attempt to repair the instrument! Please contact Miniland technical service.
6. This thermometer consists of high-quality precision parts. Do not drop the instrument. Protect it from severe impact and shock. Do not twist the instrument or the measuring sensor.
7. Please consult your doctor if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.
8. Even in the absence of fever, those who exhibit a normal temperature may still need to receive medical attention. People who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.
9. Temperature elevation may signal a serious illness, especially in adults who are old, frail, have a weakened immune system, or neonates and infants. Please seek professional advice immediately when there is a temperature elevation and if you are taking temperature for whom are:
 - Over 60 years of age (Fever may be blunted or even absent in elderly patients)
 - Having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer, chemotherapy,

- chronic steroid treatment, splenectomy)
 - Bedridden (e.g., nursing home patient, stroke, chronic illness)
 - A transplant patient (e.g., liver, heart, lung, kidney)
10. This thermometer is not intended for pre-term babies or small-for-gestational age babies. This thermometer is not intended to interpret hypothermic temperatures. Do not allow children to take their temperatures unattended.
 11. Use of this thermometer is not intended as a substitute for consultation with your physician or pediatrician.
 12. Clean the thermometer probe after each use.
 13. Do not use the thermometer on newborns or for continuous temperature monitoring purposes.
 14. Do not take a measurement while or immediately after nursing a baby.
 15. Patients should not drink, eat, or be physically active before/while taking the measurement.

4. PRODUCT INFORMATION

4.1. OVERVIEW

Infrared Thermometer measures the body temperature based on the infrared energy emitted from the eardrum or the forehead. Users can quickly get measurement results after positioning properly the temperature probe in the ear canal or forehead. Normal body temperature is a range. The following tables show that this normal range also varies by site. Therefore, readings from different sites should not be directly compared. Tell your doctor what type of thermometer you used to take your temperature and on what part of the body. Also bear this in mind if you are diagnosing yourself.

MEASURING PART	NORMAL TEMPERATURE RANGE
Forehead temperature	36.1°C - 37.5°C (97°F - 99.5°F)
Ear temperature	35.8°C - 38°C (96.4°F - 100.4°F)
Oral temperature	35.5°C - 37.5°C (95.9°F - 99.5°F)
Rectal temperature	36.6°C - 38°C (97.9°F - 100.4°F)
Axillary temperature	34.7°C - 37.3°C (94.5°F - 99.1°F)

4.2. OPERATING PRINCIPLE

In the nature, for all objects whose temperature is higher than the absolute zero (-273.15°C), for the sake of molecule's thermal motion, radiate electromagnetic wave to the surrounding ambient including the infrared wave without a break, the relationship between the density of radiation energy with the object's temperature conform to the radiation law.

The infrared thermometer's working principle is based on the fourth-power law: infer objects' radiation temperature by measuring objects' radiation infrared energy. Thermopile sensor can convert the infrared energy into thermoelectricity, and output as a detected signal after signal processing. See figure 1.

The infrared temperature sensor detects infrared energy emitted by the eardrum. A built-in lens focuses the collected energy, which is then converted into a temperature reading by the thermopiles and measurement circuits. See figure 2 and 3.

4.3. INDICATIONS FOR USE AND CONTRAINDICATIONS

The infrared thermometer takes human body temperature via the eardrum or forehead. It applies to all age groups except for babies under three months. It can be used in both home and professional environments. Do not use the thermometer if the ear is infected with otitis or suppuration.

4.4. PRODUCT STRUCTURE

4.5.1. PARTS (FIGURE 4)

The device consists of enclosure, button, a temperature sensor, Infrared temperature-measuring element, microcomputer-controlled circuit, LCD, backlight and buzzer.

- A. LCD display screen
- B. Mode/Memory button
- C. Measure/Mute-unmute button
- D. Probe
- E. Battery cover
- F. Probe cover (take it off when measuring ear temperature)

4.5.2. LCD DISPLAY (FIGURE 5)

- 1. Adult forehead mode
- 2. Child forehead mode (under 12 years old)
- 3. Ear mode
- 4. Adjust mode
- 5. Fahrenheit / Celsius degrees
- 6. Previous temperature value
- 7. Low battery indication
- 8. Mute /Un-mute icon
- 9. Current temperature value

5. INSTRUCTIONS FOR USE

5.1. FEATURES

- Magnetic probe cover and automatic mode switch between forehead & ear mode
- Quick measurement, less than 1 second
- Small Body with super large font
- Easy operation, one button design, to measure ear, forehead
- 40 sets of memories, easy to recall
- Switching between mute and un-mute mode
- Audio and visual fever warning
- Switching between °C and °F
- Auto shut-down and power-saving
- Backlight color temperature indication
- Patent technology to realize separate measurement of adult and child

5.2. BEFORE USE

5.2.1. RECOMMENDATIONS

In order to avoid the inaccuracy:

- Please make sure that the device will be used in the room only, and there is no strong conversation of wind.
- If the device is transferred from one condition to another, which has different ambient temperature, it is suggested to deposit for more than 30 minutes.
- If the tester is transferred from one condition to another, which has different ambient temperature, it is suggested to have a rest for more than 10 minutes.
- Do not hold the device for a long time, as it is highly sensitive to heat.

- The device has undergone clinical test, it is safe and accurate when using in accordance with operation manual.

5.2.2. TEMPERATURE TAKING TIPS

1. It is important to know each individual's normal temperature when they are well. This is the only way to accurately diagnose a fever. Record readings twice a day (early morning and late afternoon). Take the average of the two temperatures to calculate normal equivalent temperature. Always take the temperature in the same location, since the temperature readings may vary from different locations on the forehead.
2. A child's normal temperature can be as high as 37.7°C (99.9°F) or as low as 36.1°C (97.0°F). Please note that this unit reads 0.5°C (0.9°F) lower than a rectal digital thermometer.
3. External factors may influence ear temperatures, including when an individual has:
 - been lying on one ear or the other
 - had their ears covered
 - been exposed to very hot or very cold temperatures
 - been recently swimming or bathing

In these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature.

Use the untreated ear if prescription ear drops or other ear medications have been placed in the ear canal.

4. Holding the thermometer for too long in the hand before taking a measurement can cause the device to warm up. This means the measurement could be incorrect.
5. Patients and the thermometer should stay in steady-state room condition for at least 30 minutes.
6. Before placing the thermometer sensor onto the forehead, remove dirt, hair, or sweat from the forehead area. Wait 10 minutes after cleaning before taking measurement.
7. Use an alcohol swab to carefully clean the sensor and wait for 5 minutes before taking a measurement on another patient. Wiping the forehead with a warm or cool cloth may impact your reading. It is advised to wait 10 minutes before taking a reading.
8. In the following situations it is recommended that 3-5 measurements in the same location be taken and the highest one taken as the reading:
 - Newborn infants in the first 100 days.
 - Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.

5.2.3. INSTALLATION/REMOVAL OF BATTERIES

1. Put two AAA batteries correctly into the battery compartment.
2. Remove the batteries if the thermometer will not be used for more than two months.

Note: All memory records will be lost when uninstalling or reinstalling the battery.

5.3. TEMPERATURE TAKING INSTRUCTIONS

5.3.1. FOREHEAD TEMPERATURE MEASUREMENT (FIGURE 6)

For the forehead temperature measurement follow the next steps:

1. Press and release the Measure/Mute-unmute button in one second to power on.
2. Put the magnetic probe cover on the thermometer.
3. Press the Mode/ Memory button to choose adult mode or child mode  (for children under

12 years old)

4. Point the temperature probe at the center of the forehead, with a distance of 0-3 cm (0-1.18 inches).
5. Press and release the Measure/Mute-unmute button for 1 second. When the beep is heard, the measurement is completed.
6. Remove the thermometer from the forehead and observe the temperature.
7. The unit will shut down automatically after 13 seconds of no use.

Notes:

- The forehead measurement is an indicative reading. The measured forehead temperature can fluctuate up to 1 °F/0.5 °C from your actual body temperature. Please be aware of the factors that influence the accuracy as described in section 5.2.2, "Temperature taking tips" and section 3, "WARNINGS AND PRECAUTIONS".
- After each measurement, you can enter the recall mode and query earlier temperature readings.

WARNINGS:

- If the eyebrow area is covered with hair, sweat or dirt, please clean the area beforehand to improve the reading accuracy.
- You should wait at least 10 seconds between each measurement.
- Always check if the lens is clean.
- After each measurement, clean the temperature probe with a soft cloth, and put the thermometer in a dry and well-ventilated place.
- Always make sure the user and the thermometer have been in the same room for at least 30 minutes prior to the measurement.
- It is dangerous to make a self-diagnosis or self-treatment based on the obtained measurement results. For such purposes, please consult a doctor.

5.3.2. EAR TEMPERATURE MEASUREMENT (FIGURE 7)

For the ear temperature measurement follow the next instructions:

1. Press the Measure/Mute-unmute button to power on.
2. Remove the magnetic probe cover so the thermometer switches automatically to ear mode measurement.
 - Fit the probe snugly into the ear canal. Keep in mind that:
 - For children under 1 year: Pull the ear straight back.
3. For children aged 1 year to adult: Pull the ear up and back.
4. Press and release the Measure/Mute-unmute button in 1 second. When the beep is heard, the measurement is completed.
5. Remove the thermometer from the ear and observe the temperature.
6. The unit will shut down automatically after 13 seconds of no use.

Notes:

- Please be aware of the factors that influence the accuracy as described in section 5.2.2, "Temperature taking tips" and section 3, "WARNINGS AND PRECAUTIONS".
- After each measurement, you can enter the recall mode and query earlier temperature readings.

WARNINGS:

- Please make sure that the ear is clean, with no earwax or obstructions.
- The right ear reading may differ from the reading taken at the left ear. Therefore, always take the temperature in the same ear.
- Do not force the thermometer into the ear canal. Otherwise, the ear canal may get injured.
- When taking the temperature on an adult, gently pull the ear up and back to make sure the ear canal is straight, so that the temperature probe can receive an infrared ray from the eardrum.
- Be careful when taking the temperature on a child, whose ear canal is small.

- If the patient has an ear infection or discharge, do not use this mode.
- You should wait at least 10 seconds between each measurement.
- After each measurement, clean the temperature probe with a soft cloth, and put the thermometer in a dry and well-ventilated place.
- It is dangerous to make a self-diagnosis or self-treatment based on the obtained measurement results. For such purposes, please consult a doctor.

5.3.3. DISPLAYING THE RESULTS OF THE MEASUREMENT

- a. In forehead or ear mode.

If $32^{\circ}\text{C} \leq T \leq 37.5^{\circ}\text{C}$ ($89.6^{\circ}\text{F} \leq T \leq 99.5^{\circ}\text{F}$), the green light will be displayed, with one long beep. If $37.6^{\circ}\text{C} \leq T \leq 38.0^{\circ}\text{C}$ ($99.7^{\circ}\text{F} \leq T \leq 100.4^{\circ}\text{F}$), the orange light will be displayed with 3 short beeps, which is a warning that you may have a low fever.

If $38.1^{\circ}\text{C} \leq T \leq 42.9^{\circ}\text{C}$ ($100.6^{\circ}\text{F} \leq T \leq 109.2^{\circ}\text{F}$), the red light will be displayed with 3 short beeps, which is a warning that you may have a high fever.

- b. In adjust mode.

If $0^{\circ}\text{C} \leq T \leq 100^{\circ}\text{C}$ ($32.0^{\circ}\text{F} \leq T \leq 199.0^{\circ}\text{F}$), the white light will be displayed with one long beep.

5.4. OTHER FEATURES

5.4.1. MUTE MODE

This thermometer features an optional silent mode, allowing you to disable the beep sound during temperature measurement. To activate or deactivate this function press the Measure/Mute-unmute for 5 seconds.

When silent mode is enabled, the thermometer will operate without emitting any sound, making it ideal for use while your baby is sleeping or in quiet environments.

5.4.2. MEMORY CHECKING

There are a total of 40 memory data for recording ear and forehead temperature measurements.

To access the memory, when the thermometer is turned on or off, keep pressing the Mode/ Memory button for 5 seconds, press this button again to check the 40 sets of memories one by one. If no value, it will display "--- "

The current measurement is always stored in the last storage space. When all storage spaces have been occupied, the oldest measurement is always deleted from memory.

5.4.3. FAHRENHEIT/CELSIUS SETUP

This thermometer allows the user to switch between Celsius ($^{\circ}\text{C}$) and Fahrenheit ($^{\circ}\text{F}$). To do it, open the battery cover and use the toggle switch to change the $^{\circ}\text{C}/^{\circ}\text{F}$.

5.5. TROUBLESHOOTING

SYMPTOM	POSSIBLE CAUSE	DESCRIPTION AND SOLUTION
Failed to power on	The battery level is too low.	Replace with a new battery.
	Polarities of the batteries are reversed.	Ensure the batteries are in the right position.
	The thermometer is damaged.	Contact the MINILAND technical service.
The reading is too high	You have just come from a hot environment.	Stay in an adequately cool room for at least 30 minutes before taking a reading.
	The ambient temperature is not in range.	One long beep, and a red backlit for 3 seconds. Take a measurement under an ambient temperature between 0°C (32.0°F) and 40°C (104°F).
	Memory error or calibration is not finished. Item will be turned off automatically after 3 seconds.	One long beep, and red backlit for 3 seconds. Contact the MINILAND technical service.
	In ear or forehead mode, T>42.9°C (109.2°F). In adjust mode, T> 100°C (199°F).	One long beep, the screen will be displayed in red.
	In ear or forehead mode, T<32°C (89.6°F). In adjust mode, T<0°C (32°F).	One long beep, the screen will be displayed in red.
 The battery icon is flashing.	The power voltage is between 2.4V-2.7V.	The battery icon will flicker, it means the battery is low, but you still can use the thermometer.
 The battery icon is blinking, and the thermometer turns off automatically.	The power voltage is lower than 2.4V.	The battery icon will flicker, and the item will be turned off automatically after 3seconds. Please replace with a new battery.

6. PRODUCT CARE

Use an alcohol swab or cotton swab moistened with 70% alcohol to clean the thermometer casing and the measuring probe. After the alcohol has completely dried out, you can take a new measurement.

Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. Take care not to scratch the surface of the LCD screen.

7. TECHNICAL SPECIFICATIONS

Product name	Thermosense plus: Forehead & Ear in frared thermometer
Reference	89680
Model	AQJ-20C
Applicable regulations and laws	ISO 80601-2-56
Display	Segment LCD, 4 color LED backlight (white, green, orange, red)
Temperature units	°C / °F, switchable
Power supply	DC 3V, AAAX2
Measurement range	Measurement: 32.0°C - 42.9°C (89.6°F - 109.2°F)
	Direct mode: 0°C - 100°C (32°F - 199°F)
Accuracy (Laboratory)	Measurement mode: $\pm 0.2^{\circ}\text{C}$ / $\pm 0.4^{\circ}\text{F}$
	Direct mode: $\pm 1.0^{\circ}\text{C}$ / $\pm 2.0^{\circ}\text{F}$
Display resolution	0.1°C / °F
Automatic	13s $\pm$ 1s
Shutdown memory	40 groups of measured temperature
Operational conditions	Temperature: 5 °C-40°C(41°F-104°F)
	Humidity: 15-95%RH, non-condensing
	Atmospheric pressure:70-106kPa
Storage condition	Ambient Temperature: -20°C-55°C (-4°F-131°F)
	Relative Humidity: 0-95% RH, non-condensing
	Atmospheric pressure: 50kPa to 106kPa
Battery	2*AAA, can be used for more than 3000 times
Weight & Dimension	60g (without battery),143x35x41mm
Date of last revision of the manual	February 2025

NORMALIZED SYMBOLS

SYMBOL	DESCRIPTION
	Type BF applied part.
	Information about the manufacturer.
	Please read the instructions carefully.
	Waste electrical materials should be sent to a dedicated collection point for recycling.
	Production Date.
LOT	Batch number.
	IMPORTANT Inaccurate reading or thermometer damage may occur if the thermometer is not correctly used.
IP22	2 Protected against solid foreign objects of 12,5 mm Ø and greater; 2 If keep the thermometer in 15 degree angle, it still can prevent the water drop.
	CE marking.
	Authorized Representative In The European Community.
	Keep away from sunlight.
	Fragile, handle with care.
	Keep dry.
	This way up.

ELECTROMAGNETIC COMPATIBILITY (EMC)

EMC information-Guidance and manufacturer's declarations

CAUTION:

- The Infrared Thermometer needs special precautions regarding EMC and needs to be installed and put

ENGLISH

- into service according to the EMC information provided for in the ACCOMPANYING DOCUMENTS.
- Portable and mobile RF communications equipment can affect Infrared Thermometer.
- The Infrared Thermometer should not be used adjacent to or stacked with other equipment. Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

ELECTROMAGNETIC EMISSIONS

GUIDANCE AND MANUFACTURER'S STATEMENT - ELECTROMAGNETIC EMISSION	
The Infrared Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer should assure that it is used in such an environment.	
Emission Test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable

ELECTROMAGNETIC IMMUNITY

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY		
The Infrared Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer should assure that it is used in such an environment.		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	Not applicable	Not applicable
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM a 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM a 1 kHz
NOTE: UT is the a.c. mains voltage prior to application of the test level.		

ELECTROMAGNETIC IMMUNITY

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY				
The Infrared Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer should assure that it is used in such an environment.				
	Test Frequency (MHz)	Band (MHz)	Service	Modulation
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to wireless RF communication equipment)	385	380-390	TETRA 400	Pulse modulation 18 Hz
	400	430-470	GMRS 460 / FRS 460	FM $\pm$ 5 kHz deviation 1 kHz sine
	710	704-787	LTE Band 13,17	Pulse modulation 217 Hz
	745			
	780			
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz
	870			
	930			
	810	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1,3,4,25, UMTS	Pulse modulation 217 Hz
	870			
	930			
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450	Pulse modulation 217 Hz
			LTE Band 7	
	810	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz
	870			
	930			

ELECTROMAGNETIC IMMUNITY

GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC IMMUNITY				
The Infrared Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Thermometer should assure that it is used in such an environment.				
	Maximum Power(W)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to wireless RF communication equipment)	1.8	0.3	27	27
	2	0.3	28	28
	0.2	0.3	9	9
	2	0.3	28	28
	2	0.3	28	28
	2	0.3	28	28
	0.2	0.3	9	9

8. INFORMATION ON THE DISPOSAL OF THE BATTERIES AND THE PRODUCT

- Remember to remove the batteries and dispose of them in a responsible manner when disposing of the baby monitor.
- You should not mix products bearing the cross-out rubbish bin symbol with your domestic rubbish.
- These products should be taken to the collection points allocated by the local authorities for appropriate treatment.
- Alternatively, you may contact the point of sale where you purchased the product.

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© Miniland, S. A. 2025
Shenzhen AOJ Medical Technology Co., Ltd.
Share Info GmbH