miniland

thermotalk plus

Español	3	Português	38	Deutsch	72	Polski	106
English	21	Français	55	Italiano	89	Русский	123

ÍNDICE

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INDEX

1. Introduction
2. Contents
3. Safety instructions
4. Product features
5. Instructions for use
6. Normalized symbols
7. Battery & product disposal information

1. INTRODUCTION

Congratulations on your purchase of this Ear and forehead thermometer with pulse reading function.

The product is intended for the intermittent monitoring of human body temperature on people of all ages. This device is an ear canal IR thermometer, which means that the measure site is ear canal. The forehead temperature is obtained based on the site offset between the ear and the forehead, and the forehead temperature is only used as a reference. This thermometer is intended for home use by customer only. The patient may be an operator and can change the batteries under normal circumstance and maintain the device and its accessories according to the user manual.

Before using this device, please read all instructions in this manual carefully. For the safe and correct use of this product, please read the safety precautions explained in this user manual. Please store for future reference as it contains important information.

The features described in this user's guide are subject to modifications without prior notice.

2. CONTENTS

- 1 infrared thermometer
- 1 support
- 1 storage and travel bag
- 2 x AAA batteries
- 1 Quick guide & Guarantee card

If any of the above is missing, please contact your retailer.

3. SAFETY INSTRUCTIONS

3.1. SAFETY PRECAUTIONS

Warning marks and symbols are included for safety and a correct use of this product and prevention of any injury.

See below for the description of warning marks and symbols:

• **WARNING:** Means the possibility of personal injury in case of improper use.

• **NOTICE:** Means the possibility of injuries or damage* in the case of improper use.
[*This covers any damage to your house, personal belongings and pets]

• **PROHIBITION:** Means that the detailed items expressed in words or figures within or beside the symbol are totally forbidden. The symbol on the left means forbidden in general.

• **MUST BE OBSERVED:** This means obligatory with the detailed items expressed in words or figures within or beside the mark. The one on the left means that it is compulsory in general.

• **THE MEANING OF THE SYMBOL:** Means Type-BF device

• **The patient may be an operator and can change the batteries under normal circumstance and maintain the device and its accessories according to the user manual.**

• **It is dangerous for patients to perform a self-evaluation and self-treatment based on the results of the reading. Be sure to follow your doctor's instructions.**

**A self-evaluation may cause deterioration in the state of the patient's health.*

• **Don't touch or blow on the infra-red sensor.**

**A dirty infra-red sensor may cause an inaccurate reading.*

• **Clean the infra-red sensor with a soft wet cloth.**

** Do not clean with toilet paper or a tissue as this may scratch the infrared sensor and cause an inaccurate reading.*

• **Clean the infra-red sensor every 4 uses (following the instructions given on the manual, section 6.3. CARE AND CLEANING).**

• **Keep the machine out of children's reach**

** If a child takes his own reading, it may cause damage to his ear. Please consult the doctor at once if a battery is accidentally swallowed.*

- **When the device is used to measure the temperature of a child, it has to be operated by an adult. Adults can measure their own temperature.**
- **In the case of having performed any outdoor activities, no reading should be taken until the patient's body temperature is the same as the room temperature.**
**If the patient has used a water-filled pillow or hot water bottle, this may cause an inaccurate reading.*
- **If there is a difference in temperature between the storage area and the place where the reading is to take place, leave the thermometer at room temperature for approximately 30 minutes.**
**A failure to do so may cause an inaccurate reading.*
- **If the reading exceeds 38°C and / or the patient feels uncomfortable, the user should get in contact with a doctor.**
- **Stop using the product in the event of any pain.**
**It is likely to cause damage to the inside of the patient's ear.*
- **It is not recommended to use the product by users suffering ear disease including otitis and tympanitis.**
**The affected part is likely to get worse.*
- **Don't use the product in a wet ear.**
**It may cause damage to the ear.*
- **Avoid mechanical shock. otherwise the thermometer will be damaged.**

- **Do not dispose the batteries by throwing them into a fire.**
**They may explode.*

IMPORTANT:

- If this advice is not followed, the product guarantee may become invalid.
- The warranty excludes the user-caused failures listed below:
 - Unauthorized disassembly and modification of the product
 - If the product is dropped during application or transportation
 - An incorrect use of the product.

3.2. RECOMMENDATIONS

- The product should not be exposed to any chemical solvent, direct sunshine or high temperatures.
- Do not drop, step on or subject the product to any vibration or impact.
- Do not use a mobile phone near the thermometer when it is being used.
- For the disposal of batteries and waste, please act according to local law.
- Remove the batteries if the unit is not going to be used for a long period of time.

- Do not repair or do maintenance operations to the thermometer while it is in use.

IMPORTANT: You are advised that changes and modifications not carried out on the equipment by the MINILAND technical service are not covered by the product guarantee.

4. PRODUCT FEATURES


Intended for use:

The device is intended for the occasional reading and monitoring of human body temperature in a domestic environment.

The reference body site is the core and the mode of operation is adjusted mode.


4.1. PRODUCT STRUCTURE

4.1.1. PARTS DESCRIPTION


- | | |
|-------------------|---------------------|
| 1. Ear key | 6. Broadcast key |
| 2. Probe | 7. Forehead probe |
| 3. Support | 8. Forehead key |
| 4. LCD | 9. Battery |
| 5. Heart rate key | 10. Infrared sensor |

4.1.2. LCD DISPLAY

1. AM/PM
2. Time, temperature or pulse display
3. Memory number
4. Memory number or year display
5. Month display
6. Date display
7. Heart rate reading symbol
8. Ear temperature reading symbol
9. Speaker
10. Forehead temperature reading symbol
11. Set display
12. Body temperature reading symbol
13. Low battery symbol
14. Pulse
15. Cleaning indicator

4.2. BASIC FUNCTIONS

1. Environment temperature, current date and time display:
Range of date display falls within Gregorian calendar year 2001~2099, time display mode is 24 hours shift. LCD can display date, time and environment temperature, LCD displays time and temperature every 5 seconds.

2. Reading and display of temperature:

Readings range from 32.0°C to 42.9°C (89.6°F to 109.3°F).

3. 15 sets of readings:

Maximum of 15 sets of temperature records stored.

Temperature is stored together with date, memory number and mode (Ear or Forehead).

4. Reading display of room temperature :

Reading ranges from 5.0°C to 59.9°C (41.0°F to 139.8°F).

5. Voice prompt: readings, and/or room temperature, and/or current time.

4.3. ADDITIONAL FUNCTION

1. Function: Reading and display of heart rate.

2. Range: 40~140

3. To monitor your heart rate:

If you want to monitor your heart rate, please press [Heart rate] button, the device makes a “DI” sound, the LCD will display as in figure 2 and the device starts to read your heart rate. Meanwhile place your finger on the [Heart rate] key as in figure 1. When the device detects the heart rate, the symbol [♥] will flash to the beat of your heart.

When the reading is complete (about 15 seconds), you will hear a “DI” sound and the LCD will display your heart rate as in figure 3.

Caution:

1. Do not keep the [Heart rate] button pressed down when the device is taking a reading of your heart rate.

2. If the symbol [♥] does not flash, please make sure your finger is placed correctly and you are pressing down hard enough. Excessive pressure will

obstruct the local blood circulation and the device will not be able to detect the heart rate.

3. The maximum reading time is 15 seconds.

4. Avoid swimming or intense physical exercise at least half an hour before measuring heart rate.

5. Do not measure heart rate if your hand is shaking.

6. You must put your finger on the heart rate button (5) without moving it.

7. Whilst measuring, do not talk or move.

8. If the measurement is incorrect, you may take it again. However, you must wait 5 to 15 minutes before taking the next measurement.

9. Make sure you use the same finger if you decide to repeat the heart rate measurement so that you can compare the results.

10. The results obtained after measuring heart rate must only be regarded as a reference. They cannot be used as the basis for carrying out clinical diagnoses or evaluations. It is dangerous to draw conclusions or start treatments based on these results. Please consult your doctor and follow his indications if you don't feel well.

Figure 1


Figure 2


Figure 3

5. INSTRUCTIONS FOR USE

5.1. PREPARATIONS

1. Replace any low-voltage battery to ensure full power supply.

2. For a correct reading, the thermometer must be in the same environment 30 minutes before the reading takes place.

3. Unforeseen fluctuations in ambient temperature are likely to affect the accuracy of the measurement. When the temperature of the measurement site is different from that of the area where the thermometer was located, or when attempting to take the temperature, brake an air conditioner, it is impossible to obtain reliable results.

4. The reading must be taken in a clean ear to assure maximum accuracy.

5. The reading must be taken when the patient is in a state of rest for maximum accuracy.
6. The temperature is likely to rise immediately after exercise or bathing.

5.2. HOW TO TAKE A READING

1. To take an ear temperature:

Pull out the forehead probe from the unit, and insert the temperature probe along the eardrum as deep as possible into the ear. Press [Ⓞ] key to take the ear temperature. Upon completion you will hear two short beep sounds, the LCD will display “ooo” and the result of the reading, and you will hear the result via the speaker if the TALK function is enabled. When the “ooo” disappears from the screen, the next reading can be taken. Hold the device correctly and do not press any other keys during the ear temperature measurement.

2. To take the forehead temperature:

Put on the forehead probe onto the unit, then locate the probe gently in the center of the forehead, midway between the eyebrow and the hairline. Press [↵] key once to start the forehead measurement (the measuring time is approximately 1~2 seconds). Upon completion you will hear two short beep sounds, the LCD will display the result of the reading, and you will hear the result via the speaker if the TALK function is enabled.

3. The result will automatically be recorded in the memory. Once 15 recordings have been stored, on taking the 16th recording, the first set of data will be removed and the new recording will take the top position.
4. Perform a maximum of 3 consecutive measurements. Body temperature may not give a true reading even though the patient has recovered. For continuous measurement up to 4 times or more, you must wait at least 10 minutes before taking a new measurement.
5. Once the reading has finished, the unit will automatically return to time mode after 30 seconds.

6. The display will automatically turn off after 3 minutes after the last key has been pressed (if sleep mode is enabled).

5.3. DISPLAYING THE RESULTS OF THE MEASUREMENT

Normal body temperature is a range. A person's normal temperature range tends to decrease with age. The following table shows normal temperature ranges by age:

0 - 2 years	36.4 - 38.0 °C	97.5 - 100.4 °F
3 - 10 years	36.1 - 37.8 °C	97.0 - 100.0 °F
11 - 65 years	35.9 - 37.6 °C	96.6 - 99.7 °F
>65 years	35.8 - 37.5 °C	96.4 - 99.5 °F

The range of normal temperature varies from person to person and can be influenced by many factors such as time of day, measure from different sites, level of activity, medications... So we recommend that you practice with the thermometer on yourself and family members when you are healthy. This way you will know how the thermometer works and will feel more confident of the measurements you take when a family member is ill.


On this thermometer, the results will be showed as follows:

- A. An accurate value display for result ranges between 32.0°C and 42.9°C (89.6°F and 109.3°F)
- B. "HI" is displayed for result above 42.9°C/109.3°F
- C. "LO" is displayed for result below 32.0°C/89.6°F

5.4. ADVANCED OPERATION

5.4.1. Clock mode

1. Display time in a sequence of year, month, date, hour, minute with a flickering " : " stands for seconds. The LCD displays the current time and room temperature every 5 seconds.
 2. Press [] to show the current time and/or room temperature (except for French and Russian).
 3. Press hold [] for 3 seconds to select between english, german, french, spanish, italian, russian or no voice.
- When the unit is on voice mode, you will always hear a beep during the measurement.

4. When the TALK function is enabled, the [🔊] icon is displayed. If it is disabled, the [🔊] icon will be off. The default TALK function is enabled.

5. Press [↵] and [📺] key to retrieve the recorded temperature; press the [↵] and [📺] key for approximately 3 seconds to enter the set up mode, the time, date and LCD state (on or off) can be reset.

6. Press the [°C] and [📺] keys to select Celsius or Fahrenheit.

7. Sleep mode – The system enters sleep mode 3 minutes after the last key has been pressed.

5.4.2. SETUP MODE

Press the [↵] and [📺] keys for approximately 3 seconds to enter the setup mode.

1. Press [↵] to select between a setting of HOUR, MINUTE, YEAR, MONTH, DAY and SLEEP enable/ disable and to increase the digit.
2. Press [📺] to confirm the settings and go to the next setting.
3. Press hold [↵] to scroll through the digits quickly.

**The unit will automatically return to clock mode 30 seconds after the last key has been pressed.*

5.4.3. MEMORY MODE


1. Press [] and [] in clock mode to enter memory mode.

Press [] to display the data stored.

2. “-” will be displayed if no data is stored.

3. Maximum storage capacity is 15 sets of data with date and memory number shown.

4. The unit automatically returns to clock mode 5 seconds after the last key has been pressed.

Memory no data

Memory data

Note: Press any key to stop the sound during the announcement.

5.5. RECOMMENDATIONS

1. When using for the first time, we recommend taking the reading 3 times in the same ear. The maximum reading should be recorded in the case of different results as there may be a slight difference of between (+/- 0.3°C, +/- 0.5°F) if the thermometer is used continuously.

2. Wait only 30 seconds before taking another reading on the same patient again to avoid excessive cooling of the skin.

3. For users in sound health, both ears are almost the same in terms of body temperature. You should select the ear that shows a stable and relatively higher temperature.

4. Please clean the ear before taking a reading.

NOTE: The unit may not function correctly in an environment with electrostatic discharge so it may be necessary to reset the unit.


5.6. TROUBLES AND TROUBLESHOOTING

PROBLEMS	CHECK LIST	SOLUTIONS
No response/ Automatic reset when the insulator is taken out	Battery used up	Insert a new battery
	Battery in wrong polarity	Take out battery, reinsert battery correctly
	Poor battery contact	
Battery symbol on LCD	Low battery life	Insert a new battery
“Lo” on screen (result below 32.0°C or 89.6°F)	Thermometer placed correctly in the eardrum / forehead	Follow user’s manual to point at the eardrum / forehead.
“Hi” on screen (result below 42.9°C or 109.3°F)	Please check the operation method	Follow user’s manual for correct reading
Errp” on screen	Hardware problem	Consult with your dealer
“ErrE” on screen	Environment temperature beyond the scope of 10°C~40°C (50.0°F~104.0°F) for Ear 15°C~40°C (59.0°F~104.0°F) for Forehead	Condition the thermometer until it reaches environment temperature 10°C~40°C (50.0°F~104.0°F) for Ear 15°C~40~ (59.0°F~ 104.0°F) for Forehead
“Err” on screen	Operation error, More than 15 seconds without test result in heart rate mode	Please repeat the operation, following the instructions in user manual
Other malfunctions		Contact the MINILAND technical service.

6. MAINTENANCE

6.1. REPLACE A BATTERY

Remove the battery compartment flap. Use only new batteries (2 x AAA battery). Replace the battery compartment firmly.

6.2. MAINTENANCE, STORAGE AND CALIBRATION

Maintenance:

1. Remove any stains on the casing with a soft dry cloth.
2. Do not wash with water or detergent containing abrasive elements or benzene.

Storage:

1. Clean the thermometer and replace the cover to protect the lens after each use.
2. Do not expose the thermometer to direct sunshine, high temperatures, damp conditions, fire, flame, vibration or impact.
3. Remove the battery for long time of non use.

Calibration:

This thermometer is calibrated at the time of manufacture. If the thermometer is used according to instructions, periodic recalibration is not required. If at any time you question the accuracy of the measurement, please contact the retailer immediately.

Don't attempt to modify or reassemble the thermometer.

6.3. WARNING OF DIRT ON THE SENSOR

In order to remind the user to clean the infrared sensor, this thermometer has a special alarm. If the thermometer hasn't been used for a long time or if several measures have been taken, the LCD will display the cleaning symbol to remind the user to clean the infrared sensor. The thermometer will return back to standby mode after pressing any button.

6.4. CARE AND CLEANING

1. Keep the probe surface clean, otherwise the reading will not be accurate.
2. Clean the probe with clean paper moistened in water or alcohol, and take the reading after the water on the surface of the probe all evaporates.

Caution:

Make sure the tip of thermometer is clean. A dirty lens may cause an inaccurate reading.

7. TECHNICAL SPECIFICATION

- Product Name: Infrared thermometer
- Model: TH2071 (thermotalk plus 89068)
- Power supply: DC3V (2 x AAA battery)
- Mode of operation: Adjusted mode
- Reference body site: Core
- Temperature sensing Part: IR sensor
- Range: 32.0°C ~ 42.9°C (89.6°F ~109.3°F)
- Room Temperature: 5.0°C ~ 59.9°C (41.0°F ~139.8°F)
- Laboratory accuracy:
 - +/- 0.2°C (or +/- 0.4°F)
 - 35.5°C ~ 42.0°C (or 95.9°F ~107.6°F)
 - +/- 0.3°C (+/- 0.5°F)
 - 32.0°C ~ 35.4°C (or 89.6°F ~ 95.8°F)
 - 42.1°C ~ 42.9°C (or 107.7°F ~ 109.3°F)
- Room Temperature: +/- 2°C (4°F)
- Size: 127(L) x 39(W) x 63(D) mm
- Weight : 67 g
- Working environment: 10°C~40°C (50.0°F~104.0°F) for Ear
15°C~40°C (59.0°F ~104.0°F) for Forehead
- Relative humidity: 20%RH~85%RH
- Transportation and storage condition:
Temperature -20°C to 55°C, Relative humidity 15%-90%
- Precision/repeatability: +/- 0.3°C (+/- 0.5°F)
- Service life of the battery: With a new battery (Carbon) approximately 6 months based on a frequency of use of 5 timer per day.
- Service life of the device: 2 years
- Classification of ME equipment: Internally powered/ Continuous operation

**Dispose of used batteries responsibly, please consult your distribution center for details.*

NORMALIZED SYMBOLS

- Read the operating Instructions for use!

- Type BF Applied Part. The applied part is the probe tip.


- Disposal in accordance with Directive 2002/96/EU (WEEE)

- Complies with the European Medical Device Directive (93/42 EEC) and amended by directive 2007/47/EC requirements. Notified body TÜV Rheinland (CE0197)

- Serial number

-  • The name and the address of the manufacturer
-  • Authorized Representative in the European Community
-  • Batch code
-  • Date of manufacture

ELECTROMAGNETIC COMPATIBILITY (EMC)

- This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used
- Use of accessories other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- When the operating environment is relatively dry, strong electromagnetic interference usually occurs. At this time, the device may be affected as follows:
 - the device stops output;
 - the device turns off;
 - the device restarts;
- The above phenomenon does not affect the basic safety and essential performance of the device, and the user can use it according to the instruction. If you want to avoid the above phenomenon, please use it according to the environment specified in the manual.

Declaration - electromagnetic emission

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Emmissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device is suitable for domestic establishment and in establishment directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Immunity test	IEC 60601test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tie. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p>0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles</p>	<p>Not applicable</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>
<p>NOTE: UT is the a.c. mains voltage prior to application of the test level.</p>			

Declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>3V 0.15 MHz to 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz</p>	<p>Not applicable</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of device, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ 150kHz to 80 MHz $d = 1.2 \sqrt{P}$ 80MHz to 800 MHz $d = 2.3 \sqrt{P}$ 80MHz to 2,7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>
<p>Radiated RF IEC 61000-4-3</p>	<p>10V/m 80 MHz a 2,7 GHz</p>	<p>10V/m</p>	<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which device is used exceeds the applicable RF compliance level above, device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating device.
- b Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between
portable and mobile RF communications equipment and device**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and device, as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	0,15 MHz to 80 MHz $d = 1,2 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \cdot \sqrt{P}$	80 MHz to 2,7 GHz $d = 2,3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

8. BATTERY & PRODUCT DISPOSAL INFORMATION

- Remember to dispose of the product responsibly.
- Do not mix products bearing the crossed-out bin symbol with your general household waste.
- For the correct collection and treatment of these products take them to the collection points designated by your local authority. Alternatively contact the retailer who sold you the product.

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© Miniland, S.A. 2011
Miniland S.A. P.Ind. La Marjal I C/ La Patronal 10
Shenzhen Dongdixin Technology Co., Ltd
Shanghai International Holding Corp. GmbH (Europe)