
miniland

thermoadvanced pharma

Español	5	Português	44	Deutsch	82	Polski	122
English	25	Français	63	Italiano	103	Русский	143

Figuras · Figures


3


4

5

6

7

8

9

ÍNDICE

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INDEX

1. Introduction
2. Contents
3. Safety instructions
4. Product features
5. Instructions for use
6. Product care
7. Technical specifications
8. Information on the disposal of the batteries and the product

■ 1. INTRODUCTION

Congratulations on the purchase of this non-contact forehead digital thermometer.

This non-contact forehead thermometer is a non-sterile, reusable clinical thermometer intended for the intermittent determination of human body temperature in a touch and no touch mode on the centre of the forehead as the measurement site on people of all age. It has been carefully developed for accurate, safe, fast and contact-free human body temperature measurements (at a distance of approximately 0-5cm) in the forehead and for object temperature measurements.

This thermometer is intended for home use by customer only. The patient may be an operator and can change the batteries under normal circumstance and maintain the device and its accessories according to the user manual.

Before using this device, please read all the instructions in this manual carefully. Please keep in a safe place as it contains important information.

The features described in this user manual are subject to modifications without prior notice.

■ 2. CONTENTS

- 1 non-contact forehead thermometer
- 2x1.5V AAA batteries
- Instructions manual and warranty

If any element listed above is missing, please contact the seller.

■ 3. SAFETY INSTRUCTIONS

1. Use of this thermometer does not replace a consultation with your doctor. The measurement results can only be used as a reference. Get in touch with your doctor if you have or suspect that you may have health problems.

ENGLISH

- Users should never base their treatment or decision to seek professional help solely on the measurements made with the device, ignoring other possible illness symptoms. A delay in treatment should be avoided as it could contribute to an injury or aggravation of the underlying illness.
- Performing self-diagnoses and/or treatments based on the results of the measurements can be dangerous. Be sure to follow your doctor's instructions.
- The infrared sensor must be kept clean, dry and without deterioration at all times to ensure accurate measurements.
- Sweat or fat on the forehead can cause inaccurate measurements; make sure the forehead is clean.
- If the forehead has been cleaned, wait 5 or 10 minutes before taking the measurement.
- If the sensor is dirty, the measurement may be incorrect. Make sure it is clean before using the thermometer.
- The device is intended for occasional body temperature measurements and monitoring. Any other use is considered inappropriate and may cause material or even personal damage.
- When the device is used to measure the temperature of a child, it must be operated by an adult. Adults can measure their own temperature.
- If the temperature of the place where the thermometer is stored differs greatly from the place where the measurement is made, please wait about 30 minutes for the thermometer to equilibrate with the temperature of the room before use.
- Body temperature rises after exercising, bathing or showering, or eating. Wait at least 30 minutes before taking temperature readings.
- Do not store the thermometer at extreme temperatures, (below $-20^{\circ}\text{C}/-4^{\circ}\text{F}$ or above $55^{\circ}\text{C}/131^{\circ}\text{F}$) or extreme humidity (below 15% RH or above 90% RH). Otherwise, this may cause inaccuracies.
- Do not work with the thermometer in extreme temperatures, (below $15^{\circ}\text{C}/59^{\circ}\text{F}$ or above $40^{\circ}\text{C}/104^{\circ}\text{F}$) or extreme humidity (below 15% RH or above 85% RH). Otherwise, this may cause inaccuracies.
- Do not use the thermometer if the main casing is damaged (for example, the infrared sensor is broken). Continued use of a damaged unit may result in injury, incorrect results or serious danger.
- Please avoid mechanical shock; otherwise, the thermometer may be damaged.
- Do not attempt to open, repair or replace any part of the device except the batteries.
- Do not use any accessory that is not specified in this instruction manual.

18. Keep the thermometer out of the reach of children and pets to avoid inhalation or ingestion of small parts. Do not allow children to take their own temperature unattended.

19. Return used batteries to a suitable place for recycling. The batteries are likely to explode upon contact with fire.

20. If you are not going to use the thermometer for a long time, it is advisable to remove the batteries to avoid possible damage from chemical leaks. If the batteries are leaking, remove them carefully. Do not let the liquid touch bare skin.

21. This product needs special electromagnetic compatibility (EMC) precautions and needs to be installed and operated according to the EMC information provided. This unit may be affected by portable and mobile RF communications equipment.

22. If you have any problems with this device, contact the MINILAND technical service.

23. The patient may be the proper operator of the thermometer. Only use the thermometer as described in this user manual. Any other use is deemed improper and may result in damage to property or even personal injury. And do not service or maintain the thermometer while in use.

IMPORTANT

- Changes and modifications to the equipment not made by the MINILAND technical service are not covered by the product warranty.

- The warranty does not include the faults caused by the user indicated below:

- Faults resulting from unauthorized disassembly and modification.
- Faults resulting from accidental fall during application or transport.
- Faults resulting from not following the instructions provided in the user manual.

■ 4. PRODUCT DESCRIPTION

4.1. PARTS DESCRIPTION (Figure 1)

A. LCD screen

B. Measurement button

C. Memory button

D. Distance sensor

E. Infrared sensor

F. Battery cover

4.2. LCD SCREEN (Figure 2)

1. Forehead temperature measurement symbol

2. Object temperature measurement symbol

3. Temperature

4. Room temperature measurement symbol

5. Memory mode symbol

■ ENGLISH

6. Bluetooth symbol
7. Low battery symbol

■ 5. INSTRUCTIONS FOR USE

5.1. BEFORE USE

5.1.1. RECOMMENDATIONS

1. To obtain accurate measurements, it is very important to check the cleanliness of the infrared sensor (E) before use.

To clean the sensor (E), gently clean its surface with a cotton swab slightly moistened with alcohol and immediately dry it with a clean cotton swab. **Figure 3.**

After cleaning, let it dry for at least 20 minutes before taking temperatures.

2. Grease on the forehead can cause an inaccurate measurement, so make sure the forehead is clean. If you have had to clean the area, wait at least 5 or 10 minutes before taking the measurement.

3. If the temperature of the place where the thermometer is stored differs greatly from the place where the measurement is made, please wait about 30 minutes for the thermometer to equilibrate with the temperature of the room before use.

4. Open the battery cover and remove the insulating part when you use the thermometer for the first time; the thermometer will turn on automatically. Press any button to turn on the appliance when it is switched off. **Figure 4.**

5. Remember that in standby mode, the LCD shows the room temperature information

5.1.2. BATTERY INSTALLATION/REMOVAL

When the low battery indicator (7) appears on the LCD screen, or the thermometer stops working, you must install new batteries as soon as possible.

This product uses two 1.5 V AAA batteries; perform the following procedure to change or insert the batteries:

1. Open the battery cover. **Figure 5.**

2. Remove the old batteries and then place the new batteries in the battery compartment. Be sure to insert the batteries correctly; the positive and negative poles of the battery match the mark on the battery compartment of the thermometer. **Figure 6.**

3. Replace the battery compartment cover and press to close. **Figure 7.**

Note: To protect the environment, dispose of batteries according to local regulations and facilities. Keep batteries out of the reach of children.

5.2. FOREHEAD TEMPERATURE MEASUREMENT

External factors can influence forehead temperature, particularly when a person:

- has been exposed to very low or very high temperatures.
- has recently been swimming or bathing.
- has covered their forehead.

In these cases, the person must remain at rest and in a suitable environment for at least 30 minutes before taking the temperature.

1. Place the thermometer in your hand with your thumb on the measurement button (B). **Figure 8.**

- Do not press the button yet.
- Make sure the sensor is clean before measuring.

2. Press the measurement button (B) to start the measurement. Keep the thermometer near the forehead and in the same position until you hear two beeps. The measurement will have been completed and will be displayed on the screen. **Figure 9.**

Note: After pressing the measurement button (B), the distance sensor will detect the distance between the forehead and the thermometer.

If the distance is greater than 5 cm, you will hear a series of beeps, indicating that the thermometer is too far from the forehead. When the thermometer is brought to a correct position, the device will take the measurement, otherwise, it will return directly to standby mode.

3. Read the temperature display on the LCD screen.

To return to the main screen press the memory button (C) or wait 1 minute; the thermometer will automatically return to standby mode.

4. The thermometer will turn off automatically after 60 seconds of inactivity.

5.3. PRESENTATION OF THE RESULTS OF THE MEASUREMENT

Normal body temperature is a range. This range tends to decrease with age. The following table shows normal temperature ranges by age:

ENGLISH

0 - 2 years	38.4 - 38.0 °C	97.5 - 100.4 °F
3 - 10 years	36.1 - 37.8 °C	97.0 - 100.0 °F
11 - 65 years	35.9 - 37.6 °C	96.6 - 99.7 °F
>65 years old	35.8 - 37.5 °C	96.4 - 99.5 °F

The normal temperature range varies from person to person and can be influenced by many factors such as time of day, taking measurements at different sites, activity level, medications, etc. Therefore, we recommend that you practice using the thermometer on yourself and your family members when you are healthy. That way you will know how the thermometer works and you will feel more confident about the measurements you take when a family member is sick.

The LED indicator of this unit can be illuminated in two colours:

- Green if the measurement is between 34.0 °C (93.2 °F) and 37.9 °C (100.3 °F).
- Red if the measurement is between 38.0 °C (100.4 °F) and 42.9 °C (109.3 °F).

5.4. OBJECT TEMPERATURE READINGS

- Do not use this mode to take body temperature.
- When reading the temperature of an object, its surface temperature is displayed, which may be different from the internal temperature, especially if the surface is exposed to direct sunlight.

1. Place the thermometer in your hand and press the measurement (B) and memory (C) buttons simultaneously. The word "food" will appear on the LCD screen.
2. Place the thermometer within 5cm of the liquid or surface where you want to measure the temperature.
3. Press the measurement button (B) to start the measurement.
4. Withdraw the thermometer from the object. The LCD will show the measured temperature. You can press the memory button (C) to return to standby mode. If no button is pressed, the thermometer will automatically return to this mode after 1 minute.
5. After 1 minute of inactivity, the thermometer will turn off automatically.

5.5. MEMORY CHECK

This thermometer can memorize up to 30 forehead temperature measurements.

The current measurement is always stored in the last memory position.

When the memory is full, the oldest measurement is always deleted from the memory.

Press the memory button (C) to access the memory. You will see the storage number and temperature.

Press the memory button (C) repeatedly to access each of the measurements taken.

Press the measurement and memory buttons simultaneously to return to standby mode.

5.6. CHANGING BETWEEN FAHRENHEIT AND CELSIUS

In standby or measurement mode, you can switch between degrees Celsius (°C) and Fahrenheit (°F) by pressing the memory button (C) for a couple of seconds until you hear one or two beeps.

5.7. CONNECTION WITH EMYBABY

This thermometer connects to a smartphone or tablet via a Bluetooth connection for recording and monitoring the baby's temperatures.

To use your new infrared thermometer with your smartphone or tablet, first install the free eMyBaby application after downloading it from your Store.

When the software has been successfully installed, the eMyBaby icon will appear on the device.

Access the application and follow the steps below to record a measurement with your thermometer:

1. Make sure the Bluetooth connection of your mobile device is activated.
2. If you have not done so before, you will have to create the profile of the baby on which you want to record measurements. To do this, go to the "My Baby" section from the main menu of the eMyBaby app and follow the steps shown there to create your file.
3. Now enter the "Thermometer" section from the main menu and select the profile of the baby on which the measurement is to be carried out.
4. Select "New Measurement" and select the option to add using the "thermoadvanced pharma" thermometer.
5. The Bluetooth connection symbol (6) will appear on the thermometer screen.
6. Now measure the baby's temperature as indicated in section 5.2 of this manual.
7. After the second beep, which indicates the completion of the temperature measurement, your tablet or mobile will automatically collect the measurement taken by the thermometer.
8. Click on "Record temperature" to record the measurement and keep it in the baby's temperature history.

ENGLISH

9. You can press the memory button (C) to return to standby mode. If no button is pressed, the thermometer will automatically return to this mode after 1 minute.

10. After 1 minute of inactivity, the thermometer will turn off automatically.

5.8. PROBLEMS AND LOCATING OF FAULTS

PROBLEM	CHECKLIST	SOLUTION
No response/ Automatic reset when the insulator is removed	Batteries are dead	Replace the batteries
	Incorrect polarity of the batteries	Remove the batteries and insert them again correctly
	Poor battery contact	
The battery symbol appears on the LCD screen	Batteries low	Replace the batteries
"Lo" appears on the screen (result below 34.0 °C or 93.2 °F)	The thermometer has been positioned incorrectly	Follow the instructions of the user manual to position the thermometer correctly
"Hi" appears in the screen (result over 42.9 °C or 109.3 °F)	The operation method is incorrect	Follow the instructions of the user manual to make a correct measurement
"Err0" or "ErrE" appears on the screen	Hardware problem	Ask your distributor
"Err1" appears on the screen	The thermometer is outside the operating temperature range	Use the thermometer within the operating temperature range
	The infrared sensor temperature has not been stabilised	Wait for about 10 seconds and take the measurement again
The measurement is not accurate or there is some doubt about its result	The infrared sensor is dirty	Clean the infrared sensor as indicated in the manual

PROBLEM	CHECKLIST	SOLUTION
The measurement is not accurate or there is some doubt about its result	The measurement has been taken incorrectly	Make sure you have followed the instructions for use correctly
	The temperature of the patient and the infrared sensor have not been stabilised	Keep the thermometer and the patient in the room where the measurement is to be performed for at least 30 minutes before using it
	The measurement has been taken outdoors	Take the measurement inside a room
	The thermometer has been kept too long in the hand	Place the thermometer on a surface in the room where the measurement is to be taken and allow it to cool
The application cannot connect to the device	The Bluetooth connection has failed	Restart the app and the thermometer
Other malfunctions	Contact the MINILAND technical service.	

■ 6. CARE AND CLEANING

- To ensure the accuracy of the measurements, the distance (D) and infrared (E) sensors must be clean and free of marks. Fingerprints or other dirt may affect the accuracy of the measurements.
- Clean both sensors with a cotton swab slightly moistened with alcohol and dry them immediately with a cotton swab. To clean the thermometer casing, use a soft dry cloth.
- Never clean the thermometer with abrasive cleaning products, solvents or benzene or immerse it in water or other liquids.
- After cleaning, store it in a clean and dry place at room temperature. Never expose it to extreme temperatures and humidity or direct sunlight.

■ ENGLISH

6.1. CALIBRATION

The thermometer is calibrated during its manufacturing process. If the thermometer is used following the instructions provided in the manual, periodic re-calibration is not necessary. If at any time you have doubts about the accuracy of the measurement, immediately contact the establishment where you purchased the thermometer.

Do not attempt to modify or disassemble the thermometer.

■ 7. TECHNICAL SPECIFICATIONS

7.1. TECHNICAL SPECIFICATIONS OF THE THERMOMETER

Product	Non-contact forehead thermometer
Model	TH1009N-C (thermoadvanced pharma 89379)
Power supply	3 V DC (2 AAA batteries)
Operation mode	Adjusted mode
Body reference site	Heart
Temperature detection part	IR sensor
Body temperature measurement range	32.0°C ~ 42.9°C (89.6°F ~ 109.3°F)
Laboratory accuracy of body measurement	± 0.2 °C (or ± 0.4 °F) from 35.5 °C - 42.0 °C (or 95.9 °F - 107.6 °F) ± 0.3 °C (± 0.5 °F) for other ranges
Object temperature measurement range	-22.0°C ~ 80°C (-7.6°F ~ 176°F)
Accuracy of object measurement	± 2.0 °C (± 4.0 °F)
Room temperature measurement range	5.0°C ~ 59.9°C (41.0°F ~ 139.8°F)
Accuracy of the room measurement	± 1 °C (± 2°F)
Measuring time	1 second
Measuring interval	1 second

Weight	62 g (without batteries)
Size (length x width x depth)	130.1mm x 45.0mm x 50.3mm
Operating environment	For the forehead: 15 °C - 40 °C (59.0 °F - 104.0 °F)
	For objects: 10 °C - 40 °C (50.0 °F - 104.0 °F)
Relative humidity (RH)	15% HR ~ 85% HR
Atmospheric pressure	700hPa-1060hPa
Storage and transport conditions	Temperature -20 °C - 55 °C (-4 °F - 131 °F)
	Relative humidity: 15% - 90%
	Atmospheric pressure: 700 hPa - 1060 hPa
Maximum operating distance of the Bluetooth connection	10 m (in an open space)
Battery life	With new batteries, approximately 6 months with a frequency of use of 5 times a day
Service life	2 years
Classification of ME Equipment	Internally powered/Continuous operation
IP classification system	IP22

7.2. DESCRIPTION OF FUNCTIONS AND WIRELESS TECHNOLOGY

Frequency	2.4GHz-2.48GHz
Bluetooth® version	V4.0
Type of modulation	GFSK
No. of channels	40
Channel spacing	2MHz
EIRP	<10dBm

7.3. NORMALIZED SYMBOLS

Type BF Applied Part.

Disposal in accordance with Directive 2012/19/EU (WEEE).

Complies with the European Medical Device Directive (93/42/EEC) and amended by directive 2007/47/EC requirements.
Notified body TÜV Rheinland (CE0197)

The name and the address of the manufacturer.

Name and address of the Authorized Representative in the European Community.

Attention: Read the operating Instructions for use!

Batch code

Represents the Serial Number

Date of manufacture

IP22

The first number 2: Protected against solid foreign objects of 12,5 mm Φ and greater. The second number 2: Protected against vertically falling water drops when enclosure tilted up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15°, on either side of the vertical.

ELECTROMAGNETIC COMPATIBILITY (EMC)

- This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.
- Use of accessories other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- When the operating environment is relatively dry, strong electromagnetic interference usually occurs. At this time, the device may be affected as follows:
 - the device stops output;
 - the device turns off;
 - the device restarts;
 - the above phenomenon does not affect the basic safety and essential performance of the device, and the user can use it according to the instruction. If you want to avoid the above phenomenon, please use it according to the environment specified in the manual.

Declaration - electromagnetic emission

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for domestic establishment and in establishment directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Immunity test	IEC 60601test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to lines ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

ENGLISH

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p>0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</p> <p>0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°</p> <p>0 % UT; 250/300 cycles</p>	<p>Not applicable</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>
<p>NOTE: UT is the a.c. mains voltage prior to application of the test level.</p>			

Declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
---------------	----------------------	------------------	--

<p>Conducted RF IEC 61000-4-6</p>	<p>3V 0.15 MHz to 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz</p>	<p>Not applicable</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of device, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p>
<p>Radiated RF IEC 61000-4-3</p>	<p>10V/m 80 MHz a 2,7 GHz</p>	<p>10V/m</p>	<p>$d = 1.2\sqrt{P}$ 150KHz to 80 MHz $d = 1.2\sqrt{P}$ 80MHz to 800 MHz $d = 2.3\sqrt{P}$ 80MHz to 2,7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

ENGLISH

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which device is used exceeds the applicable RF compliance level above, device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating device.
- b Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and device, as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	0,15 MHz to 80 MHz $d = 1.2 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P}$	80 MHz to 2,7 GHz $d = 2.3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3

100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer</p> <p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

■ 8. INFORMATION ON THE DISPOSAL OF THE BATTERIES AND THE PRODUCT

- Remember to remove the batteries and dispose of them in a responsible manner when disposing of the baby monitor.
- You should not mix products bearing the crossed-out rubbish bin symbol with your domestic rubbish.
- These products should be taken to the collection points allocated by the local authorities for appropriate treatment.
- Alternatively you may contact the point of sale where you purchased the product.

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© Miniland, S.A. 2020
Miniland S.A. P.Ind. La Marjal I C/ La Patronal 10
Shenzhen Dongdixin Technology Co., Ltd
Shanghai International Holding Corp. GmbH (Europe)