miniland

eMyScale plus

eMyBaby
IS APP TO YOU!

Español	5	Português	15	Deutsch	25	Polski	35
English	10	Français	20	Italiano	30	Русский	40

eMyBaby

IS APP TO YOU!

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TABLE OF CONTENTS

1. INTRODUCTION
2. CONTENT
3. SAFETY INSTRUCTIONS
4. PRODUCT CHARACTERISTICS
5. INSTRUCTIONS FOR USE
6. MAINTENANCE
7. TECHNICAL SPECIFICATIONS
8. INFORMATION ON THE DISPOSAL OF THE BATTERIES AND THE PRODUCT

■ 1. INTRODUCTION

We would like to congratulate you on your purchase of eMyScale plus, scales for babies and children of up to 20 kg in weight, connected to your smartphone or tablet via Bluetooth using the exclusive eMyBaby application to register your baby's weight. Read the instructions manual carefully prior to using the apparatus.

Keep this manual for future use as it contains important information.

The characteristics described in this user guide are subject to change with no prior warning.

■ 2. CONTENT

1 set of scales

3 AAA batteries

1 instructions manual and guarantee

If any of the aforementioned contents are missing, contact your distributor.

■ 3. SAFETY INSTRUCTIONS

Read the safety instructions carefully prior to using this apparatus and keep it for future reference.

1. You should not strike the scales with sharp or hard objects and avoid strong vibrations.
2. Avoid heavy objects which surpass the maximum capacity of the scales (20kg/705.48oz/44lb/20000g).
3. Never wet the scales. The scales should only be cleaned with a smooth cloth and natural cleaning products.
4. Protect the scales well. If the surface is damaged the readings may be inaccurate.
5. Replace the batteries when the 'LO' indicator appears on the screen.
6. The batteries should be inserted with the poles in the correct position.
7. Remove the batteries before storing the scales for long periods of time.
8. Do not place any objects on the scale when not in use.
9. Do not place the baby on the edge of the scale to avoid any danger.
10. Protect the scale from impact with hard objects, moisture, dust, chemicals, toiletries, liquid cosmetics, temperature fluctuations and closeness to sources of heat (open fires, radiators).
11. Strong electromagnetic fields may impair the accuracy of the scale.
12. The product is for family use only.
13. Do not attempt to open or disassemble any part of the device except to replace the batteries.
14. Use the device exclusively for the purpose it is intended for.
15. Do not attempt to repair this product on your own, as opening and handling it may damage the device or expose it to hazards.

N.B: We should warn you that alterations and modifications to the equipment which are not conducted by the MINILAND technical department are not covered under the product guarantee.

■ 4. PRODUCT CHARACTERISTICS

4.1. DESCRIPTION OF THE SCALES (Figure 1)

- a. The housing
- b. Weight unit button (kg/lb/oz/g)
- c. "HOLD" button
- d. LCD screen
- e. Tare button
- f. ON/OFF button
- g. Battery compartment
- h. Brackets

4.2. DISPLAY (Figure 2)

1. Bluetooth symbol
2. "HOLD" function on/off icon
3. Minus sign
4. Battery icon
5. Tare icon
6. Digits for the weight measurement
7. Measurement unit icons

■ 5. INSTRUCTIONS FOR USE

5.1 PREPARATION

5.1.1. CHANGING THE BATTERIES

Change the batteries whenever the image on the screen is unclear and the word "LO" appears on the display.

1. Remove the cover to the battery compartment (g) located in the scales.
2. Replace the AAA type batteries. Check the batteries are correctly in place in accordance with the diagram in Figure 3.
3. Close the cover to the compartment correctly.

Note: If the scale is not going to be used for a long period of time, it is advisable to remove the batteries from the unit before storage.

5.1.2. OPERATING ENVIRONMENT

Be sure to place the scales on a firm horizontal surface in order to avoid damage to the equipment and to ensure the accuracy of the readings.

5.2 USING THE SCALES

1. Press the "ON/OFF" button (f) to turn the scales on. All the display (d) items will appear for a few seconds.
2. Press the "Weight Unit" button (b) when "0.000" appears on the screen to set the weight to "kg", "lb", "oz" or "g". The option you have selected will appear on the display (d).
3. The scale is now ready for use.
4. The Tare function enables you to remove the weight of accessories such as towels to prevent registering a false reading when weighing your baby in this manner. You should first weigh these objects and then press the Tare button (e). The word "TARE" will appear on the display and only what is placed on the scales from then on will be weighed.
5. Place the baby on the scales when "0.000" (kg/lb/oz/g) appears on the screen.
6. The definitive weight will be registered on the display when the word "HOLD" appears. The scales

ENGLISH

will turn off automatically after one minute.

7. If the baby moves during the process you may use the “HOLD” function. Turn on the scales (conducting the tare process where necessary) and place the baby in position. You will notice the reading on the display varies slightly when the baby moves. Press the “HOLD” button (c) and the word “HOLD” will start to flash on the display. As long as the word “HOLD” is not flashing the scales will read your baby’s weight correctly. The scales will turn off automatically after one minute.

Note: If you place an object on the scales before turning them on, the reading will be “0.000” when they are turned on (due to the fact the scales conduct an automatic tare when they are turned on). Remove the object and press the “Tare” button (e) to give the empty scales a reading of “0.000”.

5.3. CONNECTION WITH EMYBABY

This scale connects to a smartphone or tablet via a Bluetooth connection for recording and monitoring the baby’s weight.

Miniland has developed the eMyBaby application, available for iOS and Android, and easy to download free of charge.

In order to download the eMyBaby application you should first access the respective application store or download it directly by scanning the corresponding QR code from those illustrated below:

1. Once the application has been installed, open it. Go to LOGIN and enter your username and password to gain access, if you are already registered, or register if you have not done so before. Remember you will receive a message of confirmation via the e-mail address you provided when registering.
 2. Make sure the Bluetooth connection of your mobile device is activated.
 3. If you have not done so before, you will have to create the profile of the baby on which you want to record measurements. To do this, go to the “My Baby” section from the main menu of the eMyBaby app and follow the steps shown there to create your file.
 4. Now enter the “Scales” section from the main menu and select the profile of the baby on which the measurement is to be carried out.
 5. Select “New Measurement” and select the option to add using the “eMyScale plus” scale.
 6. The Bluetooth connection symbol (6) will appear on the scale display.
 7. Now measure the baby’s weight as indicated in section 5.2 of this manual. Your tablet or mobile phone will automatically collect the measurement taken by the scale.
 8. Click on “Record weight” to record the measurement and keep it in the baby’s weight history.
 9. After 1 minute of inactivity, the scale will turn off automatically.
- You may consult all the registers and view the progress of your baby’s weight in graph form, comparing them to WHO standards if you wish.

5.4. IDENTIFYING AND SOLVING PROBLEMS

PROBLEM	POSSIBLE CAUSES	SOLUTION
Your mobile device is unable to find eMyScale plus on taking a reading.	Your scales are not turned on.	Turn them on by pressing the "ON/OFF" button (f)
	Your mobile device is not compatible with Bluetooth BLE 4.1.	Check if your mobile device is compatible with this Bluetooth version
	Your mobile Bluetooth is not turned on.	Check the Bluetooth on your device is turned on.
	The batteries in your scales are flat.	Replace the batteries as illustrated in item 5.1.1.
	eMyScale plus is out of the Bluetooth range of your mobile device	Place your mobile device closer to the scales
The mobile device does not show the weight when weighing	You have selected the manual register option on eMyBaby	Select the scales register option
The weight reading displayed on your mobile device is incorrect	This may be due to interference during the transmission of the data	Repeat the operation
	The battery level on your scales is too low	Replace the batteries as illustrated in item 5.1.1.
The scales will not turn on	The batteries in your scales are flat.	Replace the batteries as illustrated in item 5.1.1.
	The batteries have not been put in place correctly	Insert the batteries in accordance with the instructions in item 5.1.1.
"EEE" message is showed on the display	Maximum weight of 20kg / 44lb / 705.48oz / 20000g has been exceeded	If possible remove part of the weight from the scale
"Lo" message is showed on the display	The battery level is low	Replace the batteries as illustrated in section 5.1.1.
"Unst" message is showed on the display	The scale is on an unstable surface	Place the scale on a hard, stable surface

■ 6. MAINTENANCE**6.1. CLEANING**

- Clean the surface of the scales with a smooth cloth and natural cleaning products (when necessary).
- Do not use solvents or abrasive detergents.
- Never get the scale wet or immerse it in water.

■ 7. TECHNICAL SPECIFICATIONS

- Units of weight: kg / oz / lb / g
- Maximum capacity: 20kg / 44lb / 705.48oz / 20000g

■ ENGLISH

- Minimum capacity: 0.02kg / 20g / 0.044lb / 0.7oz
- Intervals: 0.005kg / 0.01lb / 0.1oz / 5g
- Accuracy: 0.3% (of the weight of the object)
- Overload indicator: "EEEE"
- Low battery indicator: "LO"
- 1-minute automatic off function
- Bluetooth BLE 4.1.
- Batteries: 3 x AAA 1.5V

■ 8. INFORMATION ON THE DISPOSAL OF THE BATTERIES AND THE PRODUCT

- Remember to remove the batteries when disposing of this product and dispose of the same in a responsible manner.
- Do not mix the products with the crossed-out bin symbol with the rest of your domestic waste
- Once the service life of these products has ended they should be disposed of at a collection and treatment plant assigned by the local authorities. You may also contact the establishment at which you purchased the product.

Hereby, Miniland S.A. declares that the reference "89390 eMyScale plus" is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

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We declare under our own responsibility that the product:

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to which this declaration refers conforms with the relevant standards or other standardising documents:

- **Directive 2014/53/EU of the European Parliament and of the Council of 16 April 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of radio equipment and repealing Directive 1999/5/EC Text with EEA relevance.**
- ETSI EN 300 328 V2.2.2 (2019-07) Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz band; Harmonised Standard for access to radio spectrum.
- ETSI EN 301 489-17 V3.2.4:2020 (2020-09) ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems; Harmonised Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU.
- ETSI EN 301 489-1 V2.2.3 (2019-11) ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements; Harmonised Standard for ElectroMagnetic Compatibility.
- EN 62368-1:2014+A11:2017 Audio/video, information and communication technology equipment - Part 1: Safety requirements (IEC 62368-1:2014, modified).
- EN 62479:2010 Assessment of the compliance of low power electronic and electrical equipment with the basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300 GHz).
- EN IEC 61000-3-2:2019 Electromagnetic compatibility (EMC) - Part 3-2: Limits - Limits for harmonic current emissions (equipment input current \leq 16 A per phase).
- EN 61000-3-3:2013/A1:2020 Electromagnetic compatibility (EMC) - Part 3-3: Limits - Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current \leq 16 A per phase and not subject to conditional connection.
- EN 55032:2015/A11:2020 Electromagnetic compatibility of multimedia equipment - Emission Requirements.
- EN 55035:2017/A11:2020 Electromagnetic compatibility of multimedia equipment - Immunity requirements.
- **RoHS Directive (EU) 2015/863 amending Annex II to Directive 2011/65/EU on the Restriction of the use of certain hazardous substances in electrical and electronic equipment.**
- **REACH (Registration, Evaluation, Authorization and Restriction of Chemical substances) Regulation (EC) No 1907/2006.**
- **Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE).**

Miniland, S.A.
Chief Executive Officer
José Verdú Francés
Onil, Alicante, Spain.
26th November 2020

miniland
EDUCATIONAL · BABY
minilandgroup.com

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Made in China. Keep this information.

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© Miniland, S.A. 2016
Miniland S.A. P.Ind. La Marjal I C/ La Patronal 10